




MEDICAL SECTOR

OPERATION LAMPS



Kalstein®

A different accompaniment, at your service





Kalstein®

A different accompaniment, at your service

OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



KALSTEIN® Operation Lamps

Operation Lamp

Model YR02125 // YR02140



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



KALSTEIN® Operation Lamps**Operation Lamp****Model YR02125 // YR02140****Product Description:**

The quantum leap becomes especially apparent in modern operating room, with criteria for safety, efficiency and economy. Bringing expertise in surgical illumination and ceiling supply units, Kalstein offers integrated operating room systems throughout the industry for the benefit of patients and medical staff.

Model	YR02125	YR02126	YR02127	YR02128	YR02129	
Reflective treatment	Blue teak					
Dome size	70 and 60cm					
Anti reflective	Flat					
installation	Roof, wall or mobile system (only one dam for wall or mobile system)					
Dome number	1 (ceiling, wall or mobile system) 2 or 3 (ceiling)					
Horizontal oscillation	Adjustable					
Working height	Adjustable					
Work address	Adjustable					
Additional characteristics	Halogen, ceiling mount, Φ 70cm and Φ 70cm	Halogen, ceiling mount, Φ 70cm and Φ 60cm	Halogen, ceiling mount, Φ 60cm and Φ 60cm	Halogen, ceiling mount, Φ 70cm	Halogen, ceiling mount, Φ 60cm	Halogen, mobile, Φ 70cm, without battery

**#Let'sgivemore**

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



KALSTEIN® Operation Lamps

Operation Lamp

Model YR02125 // YR02140



Model	YR02130	YR02131	YR0232	YR02133	YR02136	YR02137	YR02138	YR02139	YR02140
Reflective treatment	Blue teak				Blue teak				
Dome size	70 and 60cm				70 and 60cm				
Anti reflective	Flat				Fluted				
installation	Roof, wall or mobile system (only one dam for wall or mobile system)								
Dome number	1 (ceiling, wall or mobile system) 2 or 3 (ceiling)								
Horizontal oscillation	Adjustable								
Working height	Adjustable				NIA				
Work address	Adjustable								
Additional characteristics	Halogen, mobile, Φ60cm, without battery	Halogen, mobile, Φ70cm, with battery	Halogen, mobile, Φ70cm, with battery	Halogen, mobile, Φ60cm, with battery	Halogen, ceiling mount, Φ70cm and Φ70cm	Halogen, ceiling mount, Φ70cm and Φ60cm	Halogen, ceiling mount, Φ60cm and Φ60cm	Halogen, ceiling mount, Φ70cm	Halogen, ceiling mount, Φ60cm



#Let'sgive more ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S





 **Kalstein**[®]
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



**Fundación
Jacinto Convit**

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



Humatem

Solidarité et coopération
biomédicale


<https://www.humatem.org/>



**MANIAPURE
FOUNDATION USA**

<https://www.maniapurefoundation.org/>

Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed

The 3 articles of the week!

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



1
Read this ARTICLE NOW

How does the climate affect your microscope?

Kalstein Blog

Article of the week:
How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



2
Read this ARTICLE NOW

What is the importance of the laboratory refrigerator?

Kalstein Blog

Article of the week:
What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...

23 views • 3 years ago



Kalstein's Microscopes

445 views • 1 week ago



PARIS - FRANCE
HEADQUARTER
2 Rue Jean Lantier
75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89
Cel: +33 (0) 6 80 76 07 10
sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

