




MEDICAL SECTOR

ECG DIGITAL ELECTROCARDIOGRAPH



Kalstein®

A different accompaniment, at your service





Kalstein®

A different accompaniment, at your service

OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



KALSTEIN® ECG Electrocardiograph
12 Channel Electrocardiograph Touch Screen ECG Machine
Model YR05158



#Letsgivingore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



KALSTEIN® ECG Electrocardiograph 12 Channel Electrocardiograph Touch Screen ECG Machine

Model YR05158



Product Description:

As the name suggests, the 3-channel ECG machine prints 3 waveforms / channels at a time. On 3-channel machines, the ECG signals selected by the microprocessor are amplified, filtered and sent to a 3-channel multiplexer.

Technical Description:

Model	YR05158
Power supply	Electric
Warranty	Forever
After-sales service	Return and replacement
Material	metal, plastic
Duration	1 year
Quality certification	this
Instrument Classification	Class II
Safety standard	None
Type	Cheap Price Animal ECG
Monitor	7-color LCD touchscreen Super slim design
Points of sale	Support SD card, USB disk
Analysis modes	Manual / Auto / Arrhythmia
Guides	12 Guides
Software ECGNET	Optional
Anti-baseline drift	Automatic
Interface	USB/RS232
Storage	30,000 patients save / replay
Recording	112mm *20m
Elements	specification
Characteristics	7 " display,
	touch screen
	Support Ecg analysis capabilities



#Letsgivemore ❤️

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



KALSTEIN® ECG Electrocardiograph 12 Channel Electrocardiograph Touch Screen ECG Machine

Model YR05158



Guides	12 Guides
Acquisition of potential customers	synchronously 12 leads
Input circuit	Floating; Protection circuit against the effect of the defibrillator
Input impedance	$\geq 50\text{M}\Omega$
Input circuit current	$\leq 0.05\mu\text{A}$
Recording mode	Automatic: 3CH \times 4 + 1R, 3CH \times 4, 3CHx2 + 2CHx3, 3CHx2 + 2CHx3 + 1R, 6CHx2;
	Manual: 3CH, 2CH, 3CH + 1R, 2CH + 1R;
	Rhythm: any selectable lead.
Filter	EMG Filtro: 25 Hz / 30 Hz / 40Hz/75 Hz / 100 Hz / 150Hz
	DFT Filtro: 0.05 Hz/ 0.15 Hz AC Filtro: 50 Hz / 60Hz
CMRR	>100dB
Patient current leak	<10 μA
Input circuit current	<0.05 μA
Frequent response	0.05Hz~150Hz (-3dB)
Sensitivity	2.5mm/mV, 5 mm/mV, 10 mm/mV, 20 mm/mV
Anti-baseline drift	Automatic
Constant time	$\geq 3.2\text{s}$
Noise level	<15 $\mu\text{Vp-p}$
Paper speed	12.5 mm/s, 25 mm/s, 50 mm/s
Recording mode	Thermal printing system
Paper specification	80mmx20m paper
	roll 120mmx20m paper roll



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S





 **Kalstein**[®]
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



**Fundación
Jacinto Convit**

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



Humatem

Solidarité et coopération
biomédicale


<https://www.humatem.org/>



**MANIAPURE
FOUNDATION USA**

<https://www.maniapurefoundation.org/>

Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed

The 3 articles of the week!

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



1
Read this ARTICLE NOW

How does the climate affect your microscope?

Kalstein Blog

Article of the week:
How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



2
Read this ARTICLE NOW

What is the importance of the laboratory refrigerator?

Kalstein Blog

Article of the week:
What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...

23 views • 3 years ago



Kalstein's Microscopes

445 views • 1 week ago



PARIS - FRANCE
HEADQUARTER

2 Rue Jean Lantier
75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89

Cel: +33 (0) 6 80 76 07 10

sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

