



# LABORATORY SECTOR

REFRIGERATORS AND FREEZERS



**Kalstein**®

A different accompaniment, at your service





# Kalstein®

A different accompaniment, at your service

OUR SERVICES

## Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



### Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



### Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



### #Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



### Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



### Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



### Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



**KALSTEIN® Refrigerators and Freezers**  
**368L Free Standing Deep Freezer Sprayed Steel Store Plasma Pharmacy**  
**Model YR05320**



**#Letsgivemore** ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
2 Rue Jean Lantier • 75001 Paris •  
+33 1 78 95 87 89 / +33 6 80 76 07 10 •  
<https://kalstein.eu>  
**KALSTEIN FRANCE, S. A. S**



**KALSTEIN® Refrigerators and Freezers**  
**368L Free Standing Deep Freezer Sprayed Steel Store Plasma**  
**Pharmacy**  
**Model YR05320**



**Product Description:**

Two Different Cooling System Sprayed Steel To Store Plasma And Pharmacy Standing Deep Freezer 368L. This products designed to Top refrigerator and bottom freezer, two different cooling temperature to store vaccines, blood plasma, and many other biological materials, they are widely used in researching institution and clinical sites in life science, electronic testing and medical markets.



**#Letsgivemore** ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
 2 Rue Jean Lantier • 75001 Paris •  
 +33 1 78 95 87 89 / +33 6 80 76 07 10 •  
<https://kalstein.eu>  
**KALSTEIN FRANCE, S. A. S**



**KALSTEIN® Refrigerators and Freezers****368L Free Standing Deep Freezer Sprayed Steel Store Plasma Pharmacy****Model YR05320****Technical Description:**

|                       |                            |                         |
|-----------------------|----------------------------|-------------------------|
| <b>Parameter</b>      | Model                      | YR05320                 |
| <b>Specificaitons</b> | Shelf/drawer               | Shelves/Drawer          |
|                       | Cooling method             | Direct cooling          |
|                       | Defrost                    | Manual                  |
|                       | Refrigerant                | R600a & R290            |
|                       | Power consumption(kWh/24h) | 4.15                    |
|                       | Noise(db)                  | 50                      |
|                       | Ambient temperature(°C)    | 10~32°C                 |
|                       | Temperature range(°C)      | F:-20°C~-40°C R:2°C~8°C |
| <b>Cooling system</b> | Compressor brand/QTY       | CUBIGEL/1;HUAYI/1       |
|                       | Condenser                  | Wire tube               |
|                       | Sensor                     | NTC                     |
|                       | Temperature controller     | Microprocessor          |
|                       | Display                    | Digital display         |
| <b>Power</b>          | Voltage/frequency(V/Hz)    | 220/50Hz                |
|                       | Power(W)                   | 300W                    |
|                       | Current(A)                 | 2.45A                   |
| <b>Material</b>       | Inside material            | Color sprayed steel     |
|                       | Outside material           | Color sprayed steel     |
|                       | Insulation                 | PURF                    |

**#Letsgivemore** ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
2 Rue Jean Lantier • 75001 Paris •  
+33 1 78 95 87 89 / +33 6 80 76 07 10 •  
<https://kalstein.eu>  
**KALSTEIN FRANCE, S. A. S**



**KALSTEIN® Refrigerators and Freezers**  
**368L Free Standing Deep Freezer Sprayed Steel Store Plasma**  
**Pharmacy**  
**Model YR05320**



|                    |                                    |                              |
|--------------------|------------------------------------|------------------------------|
| <b>ize</b>         | Capacity(L/cu.ft)                  | 368/12.99                    |
|                    | NT./GT.(kg)                        | 137/166                      |
|                    | Exterior size(W*D*H)(mm)           | 720*830*1870                 |
|                    | Interior size(W*D*H)(mm)           | F:580×626×625 ;              |
|                    |                                    | R:466×636×625                |
|                    | Package size(W*D*H)(mm)            | 775*920*1990                 |
|                    | CBM                                | 1.43                         |
|                    | 20GP/40GP/40HQ                     | 14/39/39                     |
| <b>Alarm</b>       | High/low temperature               | Y                            |
|                    | Power failure                      | Y                            |
|                    | Sensor error                       | Y                            |
|                    | Filter check                       | /                            |
|                    | Temperature controller failure     | Y                            |
|                    | Low battery                        | /                            |
|                    | Condenser clean                    | Y                            |
|                    | High ambient temperature           | Y                            |
|                    | Door ajar                          | Y                            |
|                    | Power failure backup system(alarm) | 8h                           |
| <b>Accessories</b> | Foot                               | Y                            |
|                    | Castor                             | Y                            |
|                    | Inside door                        | /                            |
|                    | Test hole QTY/diameter             | 2/25MM                       |
|                    | USB port                           | Y                            |
|                    | Remote alarm port                  | Optional                     |
| <b>Other</b>       | Certification                      | ISO9001,ISO13485,ISO14001,CE |



**#Letsgivemore** ❤️  
 Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
 2 Rue Jean Lantier • 75001 Paris •  
 +33 1 78 95 87 89 / +33 6 80 76 07 10 •  
<https://kalstein.eu>  
**KALSTEIN FRANCE, S. A. S**





 **Kalstein**<sup>®</sup>  
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE  
ACQUISITION  
OF A KALSTEIN  
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



**Fundación  
Jacinto Convit**

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



**Humatem**

Solidarité et coopération  
biomédicale


<https://www.humatem.org/>



**MANIAPURE  
FOUNDATION USA**

<https://www.maniapurefoundation.org/>

**Articles Kalstein**  
(Social Networks and Official Websites)



**Kalstein Blog**  
Not just be informed, be fully informed

**The 3 articles of the week!**

Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>



**1**  
Read this ARTICLE NOW

**How does the climate affect your microscope?**

**Kalstein Blog**

**Article of the week:**  
How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>



**2**  
Read this ARTICLE NOW

**What is the importance of the laboratory refrigerator?**

**Kalstein Blog**

**Article of the week:**  
What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>

**Kalstein Training Videos**  
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Alliances

23 views • 3 years ago



Kalstein's Microscopes

445 views • 1 week ago



PARIS - FRANCE  
**HEADQUARTER**

2 Rue Jean Lantier  
75001 Paris - France

**Tlf:** +33 (0) 01 78 95 87 89

**Cel:** +33 (0) 6 80 76 07 10  
sales@kalstein.eu

<https://kalstein.eu/>  
<https://www.kalstein.fr/>

