




# LABORATORY LINES

MAGNETIC STIRRER



**Kalstein**®

A different accompaniment, at your service





# Kalstein®

A different accompaniment, at your service

OUR SERVICES

## Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



### Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



### Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



### #Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



### Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



### Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



### Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



## KALSTEIN® Magnetic Stirrer

### 6-position Intelligent Magnetic Stirrer with Hotplate Model YR06140



#### Product Description

Kalstein's new lab hot plate made of Ceramic material, LED display which offers excellent chemical resistance .

- ✓ You can operate each unit separately or control all units at the same time.
- ✓ Exported Germany brushless motor for long time ,stir speed 100-2000rpm.
- ✓ LED display for Stir Speed.
- ✓ Electronic speed control for constant speed even during change in load .
- ✓ High magnetic adhesion .



<b>Model</b>	YR06140
<b>Top Plate Dimension (mm)</b>	140 x 140
<b>Top Material</b>	Black Ceramics ( square plate or round plate for each dimension anti-acid and alkali, anti-high temperature .
<b>Stir Speed</b>	100 ~ 2000rpm
<b>Stir Display</b>	Digital Display
<b>Stir speed Accuracy</b>	+/- 1 rpm.
<b>Permissible Ambient Temperature</b>	Ambient Temperature to 40°C
<b>Shell Material</b>	Superfine Aluminum alloy
<b>Voltage</b>	220v ,50Hz / 110v 60 Hz
<b>Power(w)</b>	450 x 6
<b>Motor Power for each position</b>	9W, DC14V
<b>Temperature Accuracy (°C )</b>	±1°C
<b>Stirring Flask Capacity (ml) for each position</b>	50-1000
<b>Stirring flask Capacity for water for each position</b>	50- 2000ml

Photos	Item	Quantity
	<b>Can be Removed Holder</b>	
	1. <b>Length 280mm</b> 2. <b>Diameter 10mm</b> 3. <b>Material : stainless steel .</b>	6 PCS
	<b>Clamp with screw</b>	6 PCS
	<b>Stir bar .</b>	6 PCS
	<b>User's manual</b>	1 PC



#Letsgivemore ❤️

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
2 Rue Jean Lantier • 75001 Paris •  
+33 1 78 95 87 89 / +33 6 80 76 07 10 •  
<https://kalstein.eu>  
KALSTEIN FRANCE, S. A. S





 **Kalstein**<sup>®</sup>  
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE  
ACQUISITION  
OF A KALSTEIN  
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



**Fundación  
Jacinto Convit**

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



**Humatem**

Solidarité et coopération  
biomédicale


<https://www.humatem.org/>



**MANIAPURE  
FOUNDATION USA**

<https://www.maniapurefoundation.org/>

**Articles Kalstein**  
(Social Networks and Official Websites)



**Kalstein Blog**  
Not just be informed, be fully informed

**The 3 articles of the week!**

Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>



**1**  
Read this ARTICLE NOW

**How does the climate affect your microscope?**

microscope is an instrument that allows observing and measuring objects that are too small to be seen with the naked eye. This is achieved through a complex system of lenses and mirrors, which form

**Kalstein Blog**

**Article of the week:**  
How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>



**2**  
Read this ARTICLE NOW

**What is the importance of the laboratory refrigerator?**

refrigerator is a device that is one of the most common pieces of equipment used in laboratories. It is a controlled environment (refrigerated space), various fluids and substances

**Kalstein Blog**

**Article of the week:**  
What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>

**Kalstein Training Videos**  
(Youtube Channel and Official Websites)



**Kalstein France - Jacinto Convit Foundation / Alliances**

23 views • 3 years ago



**Kalstein's Microscopes**

445 views • 1 week ago



PARIS - FRANCE  
**HEADQUARTER**  
2 Rue Jean Lantier  
75001 Paris - France

**Tlf:** +33 (0) 01 78 95 87 89  
**Cel:** +33 (0) 6 80 76 07 10  
[sales@kalstein.eu](mailto:sales@kalstein.eu)

<https://kalstein.eu/>  
<https://www.kalstein.fr/>

