



# MEDICAL SECTOR

## ANESTHETIC MACHINE



**Kalstein**®

A different accompaniment, at your service





# Kalstein®

A different accompaniment, at your service

OUR SERVICES

## Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



### Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



### Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



### #Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



### Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



### Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



### Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



**KALSTEIN® Anesthetic Machine**  
**Anesthesia Machines**  
**Model YR02064 // YR02069**



**#Letsgivemore** ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
2 Rue Jean Lantier • 75001 Paris •  
+33 1 78 95 87 89 / +33 6 80 76 07 10 •  
<https://kalstein.eu>  
KALSTEIN FRANCE, S. A. S



**KALSTEIN® Anesthetic Machine****Anesthesia Machines****Model YR02064 // YR02069****Product Description:**

Anesthesia machine has integrated all major ventilation modes such as assist, control and assist-control ventilation. It offers comprehensive monitoring and has a clear display interface, to help perform high quality anesthesia delivery and management. The anesthesia machine is suitable for inhalational anesthesia delivery for adult, pediatric and neonatal patients.

**Technical Description:**

- ✓ Anesthesia machine offers a high-precision output as low as 20ml.
- ✓ It has automatic tidal volume and compliance compensation.
- ✓ The anesthesia machine offers spontaneous and manual respiratory MV monitoring function.
- ✓ It has electronic PEEP function.
- ✓ The anesthesia delivery workstation uses autoclave and latex free components.
- ✓ Users can benefit from the tool-free maintenance and cleaning of this equipment.
- ✓ The anesthesia machine is CE compliant, so clients can feel secure in purchasing it.

**#Letsgivemore** ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
2 Rue Jean Lantier • 75001 Paris •  
+33 1 78 95 87 89 / +33 6 80 76 07 10 •  
<https://kalstein.eu>  
KALSTEIN FRANCE, S. A. S



## KALSTEIN® Anesthetic Machine

### Anesthesia Machines

Model YR02064 // YR02069



#### Technical Description:

Model	YR02064	YR02065	YR02066	YR02067	YR02068	YR02069
<b>Dimension</b>	810x860x1360 mm (WxDxH )					
<b>Weight</b>	110 kg					
<b>Type of anesthetic</b>	Halothane, enflurane, isoflurane, sevoflurane, desflurane					
<b>Gas supply</b>	O2 + N2O or O2 + N2O + Air					
<b>Inlet pressure</b>	280 ~ 650kpa (350 ~ 400kpa recommended)					
<b>Standard ventilation mode</b>	VCV, PCV, PLV, PSV, SIMV + VCV + PS, SIMV + PCV + PS					
<b>Visual display</b>	12.1 in color					
<b>Value visualization</b>	Tidal volume, minute volume, respiratory rate, I: E ratio, pressure limit					
<b>(configured and measured)</b>						
<b>Graphic display</b>	Airway pressure (Ppeak, Pmean, Pmax, Pplateau, PEEP), FiO2%,					
	Compliance (ml / cmH2O), Pressure Support Level, Inspiratory, Window, P Trigger, Inspiratory End,					
	EtCO2 and InCO2 (optional), CSI, BS%, SQI%, EMG% (optional)					
<b>Optional expansion</b>	Waveforms PT, FT, VT, PV, PF, FV loops (for compliance, resistance, leak analysis),					
	EtCO2 waveforms (when optional, EtCO2 module is present), CSI, BS%, SQI, EMG waveforms					
	(when optional CSM is present))					
<b>Additional functions</b>	12.1 "external color LCD, O2 * 2 / N2O * 2 / Air * 2, Manual / VCV / PCV / PSV / SIMV, 1x Ohmeda Tec-7 isoflurane	12.1 "external color LCD, O2 * 2 / N2O * 2 / Air * 2, Manual / VCV / PCV / PSV / SIMV, 1x Ohmeda Tec-7 isoflurane, 1xCSM	12.1 "external color LCD, O2 * 2 / N2O * 2 / Air * 2, Manual / VCV / PCV / PSV / SIMV, 1x Ohmeda Tec-7 isoflurane, 1xCSM, 1x mainstream Et CO2 (Respironics)	12.1 "external color LCD, E flowmeter, O2 * 2 / N2O * 2 / Air * 2, Manual / VCV / PCV / PSV / SIMV, 1x AGCO, 1x Ohmeda Tec-7 isoflurane	External 12.1 " color LCD, E-flowmeter, O2 * 2 / N2O * 2 / Air * 2, Manual / VCV / PCV / PSV / SIMV, 1x Ohmeda Tec-7 isoflurane, 1xCSM, 1xAGCO)	12.1 "external color LCD, E flowmeter, O2 * 2 / N2O * 2 / Air * 2, Manual / VCV / PCV / PSV / SIMV, 1x Ohmeda Tec-7 Isoflurane, 1xCSM, 1xAGCO, 1x Main Multigas Analyzer Phasein (CO2, N2O, 5 AA)



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
2 Rue Jean Lantier • 75001 Paris •  
+33 1 78 95 87 89 / +33 6 80 76 07 10 •  
<https://kalstein.eu>  
KALSTEIN FRANCE, S. A. S







 **Kalstein**<sup>®</sup>  
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE  
ACQUISITION  
OF A KALSTEIN  
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



**Fundación  
Jacinto Convit**

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



**Humatem**

Solidarité et coopération  
biomédicale

<https://www.humatem.org/>



**MANIAPURE  
FOUNDATION USA**

<https://www.maniapurefoundation.org/>



**Articles Kalstein**  
(Social Networks and Official Websites)



**Kalstein Blog**  
Not just be informed, be fully informed  
**The 3 articles of the week!**  
Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>



**1**  
Read this ARTICLE NOW  
**How does the climate affect your microscope?**  
microscope is an instrument that uses a series of lenses to magnify images of small objects. This is a very delicate instrument that requires a controlled environment to maintain its performance.  
**Kalstein Blog**  
**Article of the week:**  
How does the climate affect your microscope?  
Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>



**2**  
Read this ARTICLE NOW  
**What is the importance of the laboratory refrigerator?**  
Refrigerator is a piece of laboratory equipment that is used to maintain a controlled environment (refrigerated space), without fluctuations in temperature.  
**Kalstein Blog**  
**Article of the week:**  
What is the importance of the laboratory refrigerator?  
Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>

**Kalstein Training Videos**  
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...

23 views • 3 years ago



Kalstein's Microscopes

445 views • 1 week ago



PARIS - FRANCE  
**HEADQUARTER**

2 Rue Jean Lantier  
75001 Paris - France

**Tlf:** +33 (0) 01 78 95 87 89

**Cel:** +33 (0) 6 80 76 07 10

[sales@kalstein.eu](mailto:sales@kalstein.eu)

<https://kalstein.eu/>  
<https://www.kalstein.fr/>

