




LABORATORY LINES

BIOLOGICAL MICROSCOPE



Kalstein®

A different accompaniment, at your service





Kalstein®

A different accompaniment, at your service

OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



Biological Microscope

Model YR0239



 **Kalstein**
Binocular Microscope

YR0239



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
Optimum Business Center 450 Rue Baden Powell,
• 34000 Montpellier, France.
Tél: +33 467158849 / +33 680760710 / +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A. S






Biological Microscope

Model YR0239

Description:

- ✓ Sharp Image Display with Excellent Optical Function.
- ✓ Small Space Occupancy with Integral Stand Design, Low Environment Requirement with Anti-mould Technology.
- ✓ Comfortable Operation with Ergonomic Design Principle.



Technical specifications	
Optical System	Finite Optical System
Eyepiece	WF10×/18
Viewing Head	Monocular Head Inclined at 30°
	Seidentopf Binocular Viewing Head Inclined at 30°, Interpupillary 47-78mm
	Seidentopf Trinocular Viewing Head Inclined at 30°, Interpupillary 47-78mm
	Digital Viewing Head Inclined at 30°, Interpupillary 47-78mm, 5M pixels, Support Wifi
Achromatic Objective	4× NA=0.1 WD=18
	10× NA=0.25 WD=7
	40× NA=0.65 WD=0.53
	100× NA=1.25 WD=0.13
Nosepiece	Backward Quadruple Nosepiece
Stage	Double Layers Mechanical Stage 140X140mm, Moving Range 75X50mm
Condenser	NA1.25 Abbe Condenser
	Dark Field Condenser (Dry, Oil)
Focusing System	Coaxial Coarse and Fine Adjustment, Fine Division 0.002mm, Moving Range 25mm
Illumination	Halogen Lamp 6V/20W
	S-LED Illumination
	S-LED Illumination with Recharger
	Kohler Illumination
	ECO Function
	Plane-concave Mirror
Fluorescent Attachment	Fluorescent Illumination, LED lamp
	Fluorescent Illumination, 100W HBO lamp
Polarization Set	
Phase Contrast Kit	



#Letsgivemore


Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
Optimum Business Center 450 Rue Baden Powell,
• 34000 Montpellier, France.
Tél: +33 467158849 / +33 680760710/ +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A. S



Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed
The 3 articles of the week!
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



1
Read this ARTICLE NOW
How does the climate affect your microscope?
A microscope is an instrument that allows observing and measuring objects that are too small to be seen with the naked eye. This is achieved through a complex system of lenses, which form
Kalstein Blog
Article of the week:
How does the climate affect your microscope?
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



2
Read this ARTICLE NOW
What is the importance of the laboratory refrigerator?
A refrigerator in a laboratory is one of the most important equipment for scientists to maintain, in a controlled environment (refrigerated space), various fluids and substances.
Kalstein Blog
Article of the week:
What is the importance of the laboratory refrigerator?
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Alliances
23 views • 3 years ago



Kalstein's Microscopes
445 views • 1 week ago



MONTPELLIER - FRANCE
HEADQUARTERS
Optimum Business Center,
450 Rue Baden Powell,
34000 Montpellier - France.

Fax: +33 (0) 4 67 15 88 49
Tel: +33 (0) 6 80 76 07 10
sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

Official Websites for Latin America.

