




Kalstein®

A different accompaniment, at your service

REAGENTS & CONSUMABLES SECTOR

HISTOPATHOLOGY ACCESORIES AND
CONSUMABLES





Kalstein®

A different accompaniment, at your service

OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



KALSTEIN® Reagents & Consumables
Histopathology Accessories and Consumables



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S





Disposable Microtome Blade:

YRA18: The disposable microtome blade can use for the ultrathin section an efficient section on tissue slice at Pathology & Histology. The disposable microtome blade can work both Rotary and Sledge type of Microtome, and the blade is also easy to replace. The plastic blade box is equipped with a handle to help user to take the blade from the box. The use of blade for slicing operations need to be very careful, the cutting edge of the blade is very sharp. Size of each blade: 80mm*8mm*0.25mm

YRA19: The disposable microtome blade can use for the ultrathin section an efficient section on tissue slice at Pathology & Histology. The disposable microtome blade can work both Rotary and Sledge type of Microtome, and the blade is also easy to replace. The plastic blade box is equipped with a handle to help user to take the blade from the box. The use of blade for slicing operations need to be very careful, the cutting edge of the blade is very sharp. Size of each blade: 80mm*14mm*0.31mm



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
 2 Rue Jean Lantier • 75001 Paris •
 +33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S





YRA24 Embedding Mold:

Size: 7mm*7mm*5mm

15mm*15mm*5mm

24mm*24mm*5mm

30mm*24mm*5mm

37mm*24mm*5mm

YRA32 Forceps:
Material: Stainless steel
Length: 16cm



#Let'sgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S





Block Paraffin:

YRA33: Melting Point
 50°C-52°C, 52°C-54°C, 54°C-56°C, 56°C-58°C, 58°C-60°C, 60°C-62°C, 62°C-64°C

Microscopes Slides:

YRA37: 50pcs/box Clear glass Ground edges 25.4*76.2mm(1"*3"), 1mm-1.2mm thick



Plastic Cassette With Cover:

YRA21: Snap-latch and one-piece with attached lid allows for one-hand operation. The lids can be opened and closed as often as possible. 45°angled writing surface.



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
 2 Rue Jean Lantier • 75001 Paris •
 +33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S





 **Kalstein**[®]
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



**Fundación
Jacinto Convit**

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



Humatem

Solidarité et coopération
biomédicale


<https://www.humatem.org/>



**MANIAPURE
FOUNDATION USA**

<https://www.maniapurefoundation.org/>

Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed

The 3 articles of the week!

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



1
Read this ARTICLE NOW

How does the climate affect your microscope?

Kalstein Blog

Article of the week:
How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



2
Read this ARTICLE NOW

What is the importance of the laboratory refrigerator?

Kalstein Blog

Article of the week:
What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...

23 views • 3 years ago



Kalstein's Microscopes

445 views • 1 week ago



PARIS - FRANCE
HEADQUARTER

2 Rue Jean Lantier
75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89

Cel: +33 (0) 6 80 76 07 10

sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

