




# LABORATORY LINES

ULTRASONIC BATHS



**Kalstein**®

A different accompaniment, at your service





# Kalstein®

A different accompaniment, at your service

OUR SERVICES

## Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



### Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



### Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



### #Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



### Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



### Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



### Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.




## KALSTEIN® Ultrasonic Baths Carburator Ultrasonic Cleaner



### Product Description

- LCD display and control, simple operation;
- One-key degassing can quickly remove air in water, and the ultrasonic cleaning effect is better;
- Independently developed MCU-sweep ultrasonic generator drive, uniform and strong ultrasonic effect;
- Independently developed industrial high Q value transducer with high ultrasonic conversion efficiency and long service life;
- Power adjustment function, can meet different cleaning requirements;
- Imported high strength glue, seedless nail bonding process, higher ultrasonic conversion efficiency;
- Digital timing, cleaning time 10s-100min and digital heating control, temperature 20-80°C arbitrary setting;
- The inner groove adopts imported SUS304, 1.0mm stainless steel stamping, ultrasonic exchange effect is better, longer service life;
- The washing basket is made of high quality 304 stainless steel mesh welding, surface electrolytic polishing treatment;
- The shell is made of high quality stainless steel, with better anticorrosion effect.



LCD  
LCD display



Temperature adjustable



Degassing



Power adjustable



Heating



Timing



Deep cleaning



Industrial ultrasonic transducer

Model	Capacity	Working Frequency	Ultrasonic power	Degas	Heating Power	Timing time	Temperature range	Power Regulation	Cover	Drain Valve	Washing Basket	Inner Tank Dimension	Packing Size	Gross weight
YR09120	3.2L	39 ± 1.5KHz	120W	Yes	100W	10 100	20-80°C	0-100%	Yes	Yes	optional	240*135*100mm	342*240*355mm	4.77KG
YR09120L	5.0L		120W	Yes	200W			0-100%			optional	240*135*150mm	342*240*355mm	5.4KG
YR09180	4.5L		180W	Yes	200W			0-100%			optional	300*150*100mm	402*252*355mm	5.65KG
YR09180L	6.8L		180W	Yes	200W			0-100%			optional	300*150*150mm	402*252*355mm	6.2KG
YR09240L	11L		240W	Yes	300W			0-100%			optional	300*240*150mm	412*357*380mm	8.5KG
YR09360	15L		360W	Yes	400W			0-100%			optional	330*300*150mm	472*437*420mm	11.5KG
YR09360L	20L		360W	Yes	400W			0-100%			optional	330*300*200mm	472*437*470mm	12.5KG
YR09480	22L		480W	Yes	500W			0-100%			optional	500*300*150mm	675*454*430mm	14.5KG
YR09600	30L		600W	Yes	600W/220v 500w/110v			0-100%			optional	500*300*200mm	657*450*475mm	16.0KG

(voltage):220V±10% 50/60Hz 110V±10% 50/60Hz



#Letsgivingmore ❤️

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
2 Rue Jean Lantier • 75001 Paris •  
+33 1 78 95 87 89 / +33 6 80 76 07 10 •  
<https://kalstein.eu>  
KALSTEIN FRANCE, S. A. S







 **Kalstein**<sup>®</sup>  
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE  
ACQUISITION  
OF A KALSTEIN  
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



Fundación  
Jacinto Convit

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



**Humatem**

Solidarité et coopération  
biomédicale

<https://www.humatem.org/>




**MANIAPURE  
FOUNDATION USA**

<https://www.maniapurefoundation.org/>



**Articles Kalstein**  
(Social Networks and Official Websites)



**Kalstein Blog**  
Not just be informed, be fully informed

**The 3 articles of the week!**

Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>



**1**  
Read this ARTICLE NOW

**How does the climate affect your microscope?**

microscope is an instrument that allows observing and measuring objects at a very small scale. This is achieved through a complex system of lenses, which form

**Kalstein Blog**

**Article of the week:**  
How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>



**2**  
Read this ARTICLE NOW

**What is the importance of the laboratory refrigerator?**

refrigerator is a device that is one of the most used equipment in laboratories to maintain, in a controlled environment (refrigerated space), various fluids and substances

**Kalstein Blog**


**Article of the week:**  
What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>

**Kalstein Training Videos**  
(Youtube Channel and Official Websites)



**Kalstein France - Jacinto Convit Foundation / Allianc...**  
23 views • 3 years ago



**Kalstein's Microscopes**  
445 views • 1 week ago



PARIS - FRANCE  
**HEADQUARTER**  
2 Rue Jean Lantier  
75001 Paris - France

**Tlf:** +33 (0) 01 78 95 87 89  
**Cel:** +33 (0) 6 80 76 07 10  
[sales@kalstein.eu](mailto:sales@kalstein.eu)

<https://kalstein.eu/>  
<https://www.kalstein.fr/>

