


Kalstein A different accompaniment, at your service

LABORATORY SECTOR

REFRACTOMETER





All rights reserved ® KALSTEIN France S. A. S., 2 Rue Jean Lantier • 75001 Paris • +33 1 78 95 87 89 / +33 6 80 76 07 10 • https://kalstein.eu KALSTEIN FRANCE, S. A. S





OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.





#Letsgivemore **

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.





Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or.air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



KALSTEIN® Refractometer

Clinical Refractometer

Model YR05915







#Letsgivemore 💗

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved ® KALSTEIN France S. A. S., 2 Rue Jean Lantier • 75001 Paris • +33 1 78 95 87 89 /+33 6 80 76 07 10 • https://kalstein.eu

KALSTEIN FRANCE, S. A. S



KALSTEIN® Refractometer

Clinical Refractometer

Model YR05915





Product Description:

Designed for clinical usage and intended for both Veterinary and Human. This unit provides quick and accurate indication of vital fluids levels. It has triple scale that gives Urine Specific Gravity(SG), Total Serum Protein(SP) and Refractive Index (RI). It's FULLY EQUIPPED with ATC that automatically adjust itself to correct temperature discrepancies during use. Capable of measuring with high accuracy and provides precise results. It's equipped with vivid and sharp reticle chart for easy and comfortable reading. It uses ambient light only and no batteries are needed. Made by experienced craftsmen and are unlikely to break down!

Features:

- ✓ High quality body construction
- ✓ Uses ambient light only which means battery or power source is not required
- ✓ Easy to focus and calibrate
- ✓ Cushioned with soft & comfortable non-slip rubber
- ✓ With ATC function(Auto Temperature Compensation)
- ✓ Durable and built to last long



Technical Specifications:

Model	YR05915
Description	ATC Clinical 0- 12g/dl 1.000-1.050sg 1.3300RI-1.3600RI optical refractomete
Measurement Range	0-12g/dl 1.000-1.050sg 1.3300Rl-1.3600Rl
Min. Div.	0.2g/dl 0.001 SG 0.001 RI
Accuracy	±0.2g/dl ±0.001 SG ±0.001 RI
ATC	Yes
ATC Temp.	10-30° C
Material	Pure aluminum



#Letsgivemore 💗

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved ® KALSTEIN France S. A. S. 2 Rue Jean Lantier • 75001 Paris • +33 1 78 95 87 89 /+33 6 80 76 07 10 • https://kalstein.eu

KALSTEIN FRANCE, S. A. S







#LETSGIVEMORE >

WITH THE ADCQUISITION OF A KALSTEIN EQUIPMENT

YOU MAKE A CONTRIBUTION TO:





https://onetreeplanted.org



Solidarité et coopération biomédicale

https://www.humatem.org/



https://www.maniapurefoundation.org/



ANY DOUBT?Contact us!

Articles Kalstein

(Social Networks and Official Websites)







Article of the week: How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites https://kalstein.eu/articles





Article of the week: What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites https://kalstein.eu/articles



PARIS - FRANCE HEADQUARTER

2 Rue Jean Lantier 75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89 **Cel:** +33 (0) 6 80 76 07 10 sales@kalstein.eu

https://kalstein.eu/ https://www.kalstein.fr/

Kalstein Training Videos

(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...

23 views • 3 years ago



Kalstein's Microscopes

445 views • 1 week ago

