


Kalstein A different accompaniment, at your service

MEDICAL SECTOR

ELECTROSURGICAL UNIT





All rights reserved ® KALSTEIN France S. A. S., 2 Rue Jean Lantier • 75001 Paris • +33 1 78 95 87 89 /+33 6 80 76 07 10 • https://kalstein.eu KALSTEIN FRANCE, S. A. S





OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.





#Letsgivemore **

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.





Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or.air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



KALSTEIN® Electrosurgical Unit

Electrosurgical Generators

Model YR02145 // YR02146







#Letsgivemore 💗

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved ® KALSTEIN France S. A. S., 2 Rue Jean Lantier • 75001 Paris • +33 1 78 95 87 89 /+33 6 80 76 07 10 • https://kalstein.eu

KALSTEIN FRANCE, S. A. S



KALSTEIN® Electrosurgical Unit

Electrosurgical Generators

Model YR02145 // YR02146





Product Description:

YR electrosurgical series generator is a versatile electrosurgical unit for all around surgical procedures through all types of tissue. It has integrated seven power output, tissue responding and automatic solid or split return electrode recognition functions, with which surgeons can perform highly precise and safe operations.

Application:

YR is applicable in various surgeries such as general surgery, dermatology, plastic surgery, endoscopy, vascular surgery, heart/thoracic surgery, minimally invasive surgery, neurosurgery and ophthalmology.

Technical Data:

Model	YR02145		YR02146		
Watts	350w		150w		
Monopolar					
Cutting mode	Pure	Mix 1	Mix 2	Mix 3	
Max. Output power (watts)	350w (gray) / 150w (white)	300w (gray) / 120w (white)	200w (gray) / 80w (white)	150w (gray) / 60w (white)	
Load (ohms)	500				
Open circuit Vp-p	1,800 / 1,200v	2,700 / 2,000v	2,400 / 1,700v	2,900 / 3,000v	
Crest factor	Sun May 01 2022 00:00:00 GMT + 0200 (Central European Summer Time)	Thu Sep 01 2022 00:00:00 GMT + 0200 (Central European Summer Time)	Wed Mar 02 2022 00:00:00 GMT + 0100 (Central European Standard Time)	Mon Jan 03 2022 00:00:00 GMT + 0100 (Central European Standard Time)	
Coagulation mode					
Point	100w (gray) / 70w (white)	500	4,600 / 3,900v	Wed May 04 2022 00:00:00 GMT + 0200 (Central European Summer Time)	
Gentle	100w (gray) / 70w (white)	500	3,600 / 3,300v	Wed May 04 2022 00:00:00 GMT + 0200 (Central European Summer Time)	



#Letsgivemore 💗

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved @ KALSTEIN France S. A. S 2 Rue Jean Lantier • 75001 Paris • +33 1 78 95 87 89 / +33 6 80 76 07 10 • https://kalstein.eu

KALSTEIN FRANCE, S. A. S



KALSTEIN® Electrosurgical Unit

Electrosurgical Generators







	Connecting accessories and supplies			
Finger pencil	one			
Foot switch pencil	two			
Patient return electrode	one			
Foot switch	one			
AC power supply	one			
Potential equalization	one			
Surgical smoke evacuator	1 (but surgical fume extractor is optional)			
Power characteristic				
Rate frequency	512 kHz (110 / 240vAC, 60 / 50Hz)			
Isolated output	AND IT IS			
Quantitative visualization of the contact area	YES, for the «Split» type of neutral electrode			
Contact area indicator lamp	10 lamps (one for 10% contact), alarm beep and flash if contact <20%			
Definition of Solid / Split Neutral Electrode	YES, the 'Split' or 'Solid' indicator light comes on			
Adjustable alarm volume	YES			
Operating duty cycle	10 seconds ON, 30 seconds OFF			
cooling method	Convection			
	Regulation			
Type according to EN 60601-1	CF			
Protection class according to EN 60601-	class I			
Classification according to MDD 93/42 / EEC	IIb			
	Physical			
Dimension (WxDxH)	330 x 450 x 150 mm			
Weight	7.0 kg			



#Letsgivemore 💗

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved ® KALSTEIN France S. A. S., 2 Rue Jean Lantier • 75001 Paris • +33 1 78 95 87 89 / +33 6 80 76 07 10 • https://kalstein.eu KALSTEIN FRANCE, S. A. S







#LETSGIVEMORE >

WITH THE ADCQUISITION OF A KALSTEIN EQUIPMENT

YOU MAKE A CONTRIBUTION TO:





https://onetreeplanted.org



Solidarité et coopération biomédicale

https://www.humatem.org/



https://www.maniapurefoundation.org/



ANY DOUBT?Contact us!

Articles Kalstein

(Social Networks and Official Websites)







Article of the week: How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites https://kalstein.eu/articles





Article of the week: What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites https://kalstein.eu/articles



PARIS - FRANCE HEADQUARTER

2 Rue Jean Lantier 75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89 **Cel:** +33 (0) 6 80 76 07 10 sales@kalstein.eu

https://kalstein.eu/ https://www.kalstein.fr/

Kalstein Training Videos

(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...

23 views • 3 years ago



Kalstein's Microscopes

445 views • 1 week ago

