




# LABORATORY SECTOR

MOISTURE ANALYZER



**Kalstein**®

A different accompaniment, at your service





# Kalstein®

A different accompaniment, at your service

OUR SERVICES

## Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



### Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



### Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



### #Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



### Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



### Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.




### Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



**KALSTEIN® Moisture Analyzer**  
**Halogen Moisture Analyzer**  
Model YR05558 // YR05561



**#Letsgivemore** ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
2 Rue Jean Lantier • 75001 Paris •  
+33 1 78 95 87 89 / +33 6 80 76 07 10 •  
<https://kalstein.eu>  
KALSTEIN FRANCE, S. A. S



## KALSTEIN® Moisture Analyzer Halogen Moisture Analyzer Model YR05558 // YR05561



### Product Description:

- ✓ 7 inch touch screen
- ✓ Automatic door/Transparent window
- ✓ Highlight halogen lamp
- ✓ Polymer materials/Heat resistant molding
- ✓ Dual motor drive/Double fan design
- ✓ Temp. Control/Stainless steel chamber
- ✓ Real-time dynamic curve/Date and Time
- ✓ Temperature/Time set
- ✓ The percentage of moisture content
- ✓ Store historical set
- ✓ Connect PC and printer
- ✓ Option: Printer

### Technical Description:

Model	YR05558	YR05559	YR05560	YR05561
Capacity	110g		210g	
Readability	0.001g			
Operating Temp	5°C-35°C			
Pan Size	Ø100mm			
Heating Source	Halogen Lamp			
Temp.Sensor	PT100			
Temp.Range	40°C-230°C			
Temp.Step	1°C			
Moisture Range	0.00%-100.00%			
Moisture Readability	0.01%			
Dry Range	100.00%-0.00%			
Dry Readability	0.01%			
Temp Setting	40°C-230°C By 1°C Step			
Time Setting	1-99min By 10s			
Number Of Storage	20			
Number Of History	20			
Heater	220V ±15% 50HZ/110V±15% 60HZ/400W			
Balance	Input:220V 50HZ/110V 60HZ Output:9V			
Machine	Input: 220V 50HZ/110V 60HZ Output:24V Dual-drive			
Host Dimension	215x195x415mm			
Indicator Dimension	200x75x135mm			
Operation Height	275mm			
The Chamber Height	35mm			
Display	7 inch touch panel			
Interface	RS232/RJ45/USB...			
Communication Mode	Cable	Wireless	Cable	
Rechargeable Battery	No	10,000mAh	No	



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
2 Rue Jean Lantier • 75001 Paris •  
+33 1 78 95 87 89 / +33 6 80 76 07 10 •  
<https://kalstein.eu>  
KALSTEIN FRANCE, S. A. S







 **Kalstein**<sup>®</sup>  
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE  
ACQUISITION  
OF A KALSTEIN  
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



**Fundación  
Jacinto Convit**

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



**Humatem**

Solidarité et coopération  
biomédicale

<https://www.humatem.org/>




**MANIAPURE  
FOUNDATION USA**

<https://www.maniapurefoundation.org/>



**Articles Kalstein**  
(Social Networks and Official Websites)



**Kalstein Blog**  
Not just be informed, be fully informed

**The 3 articles of the week!**

Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>



**1**  
Read this ARTICLE NOW

**How does the climate affect your microscope?**

**Kalstein Blog**

**Article of the week:**  
How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>



**2**  
Read this ARTICLE NOW

**What is the importance of the laboratory refrigerator?**

**Kalstein Blog**

**Article of the week:**  
What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>

**Kalstein Training Videos**  
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...

23 views • 3 years ago



Kalstein's Microscopes

445 views • 1 week ago



PARIS - FRANCE  
**HEADQUARTER**

2 Rue Jean Lantier  
75001 Paris - France

**Tlf:** +33 (0) 01 78 95 87 89

**Cel:** +33 (0) 6 80 76 07 10

[sales@kalstein.eu](mailto:sales@kalstein.eu)

<https://kalstein.eu/>  
<https://www.kalstein.fr/>

