




# LABORATORY SECTOR

## INCUBATORS



**Kalstein**®

A different accompaniment, at your service





# Kalstein®

A different accompaniment, at your service

OUR SERVICES

## Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



### Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



### Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



### #Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



### Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



### Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



### Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



**KALSTEIN® Incubators**  
**Heating Incubator**  
**Model YR0205 // YR0206**



**#Letsgivemore** ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
2 Rue Jean Lantier • 75001 Paris •  
+33 1 78 95 87 89 / +33 6 80 76 07 10 •  
<https://kalstein.eu>  
**KALSTEIN FRANCE, S. A. S**



## KALSTEIN® Incubators

### Heating Incubator

Model YR0205 // YR0206



#### Product Description:

- ✓ Interior Volume: 65L, 125L
- ✓ Auxiliary Temp. control
- ✓ Leakage Protection
- ✓ Temperature Range: RT+ 5°C~66°C
- ✓ Temperature Fluctuation:  $\pm 0.5^{\circ}\text{C}$
- ✓ Temperature Uniformity:  $\pm 1^{\circ}\text{C}$
- ✓ Temperature Resolution:  $0.1^{\circ}\text{C}$

#### Options:

- ✓ Multi-segment programmable control
- ✓ Built-in printer
- ✓ RS485 interface
- ✓ UV sterilizer

#### Safety:

- ✓ Equipped with leakage protection.
- ✓ Equipped with spare temperature control which ensures the product work normally even the main temp. control failed.

#### Key Features:

- ✓ Unique air duct design, and good temperature uniformity.
- ✓ Microprocessor controller(with temperature correction and timing function).
- ✓ Large LCD screen display.
- ✓ High-quality stainless steel chamber, removable shelf, easy to clean.
- ✓ Silicon sealing ring for reliable sealing.
- ✓ Optional printer or RS485 interface which can print or connect the computer to realize remote control and alarm.
- ✓ Anti-hot handle

#### Technical Parameters:

- ✓ Temp. Fluctuation:  $\pm 0.5^{\circ}\text{C}$ .
- ✓ Temp. Uniformity:  $\pm 1^{\circ}\text{C}$
- ✓ Temp. Resolution:  $0.1^{\circ}\text{C}$
- ✓ Power Supply: AC 110V, 50Hz



Model	Temp. Range	Exterior Size (W×D×H)cm	Package Size (W×D×H)cm	Shelf	Power Rating (W)
65LTS	RT+ 5°C~66°C	52x48x66	36x34x45	2	110
125LTS		63x57x76	50x43x56	2	110



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
 2 Rue Jean Lantier • 75001 Paris •  
 +33 1 78 95 87 89 / +33 6 80 76 07 10 •  
<https://kalstein.eu>  
 KALSTEIN FRANCE, S. A. S





 **Kalstein**<sup>®</sup>  
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE  
ACQUISITION  
OF A KALSTEIN  
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



**Fundación  
Jacinto Convit**

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



**Humatem**

Solidarité et coopération  
biomédicale


<https://www.humatem.org/>



**MANIAPURE  
FOUNDATION USA**

<https://www.maniapurefoundation.org/>

**Articles Kalstein**  
(Social Networks and Official Websites)



**Kalstein Blog**  
Not just be informed, be fully informed

**The 3 articles of the week!**

Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>



**1**  
Read this ARTICLE NOW

**How does the climate affect your microscope?**

**Kalstein Blog**

**Article of the week:**  
How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>



**2**  
Read this ARTICLE NOW

**What is the importance of the laboratory refrigerator?**

**Kalstein Blog**

**Article of the week:**  
What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>

**Kalstein Training Videos**  
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...

23 views • 3 years ago



Kalstein's Microscopes

445 views • 1 week ago



PARIS - FRANCE  
**HEADQUARTER**

2 Rue Jean Lantier  
75001 Paris - France

**Tlf:** +33 (0) 01 78 95 87 89

**Cel:** +33 (0) 6 80 76 07 10

[sales@kalstein.eu](mailto:sales@kalstein.eu)

<https://kalstein.eu/>  
<https://www.kalstein.fr/>

