




LABORATORY LINES

HEATING MANTLES



Kalstein®

A different accompaniment, at your service





Kalstein®

A different accompaniment, at your service

OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



KALSTEIN® Heating Mantles

Series YR05021



Product Description

YR heating mantles are the most widely used heating mantles in laboratories. The maximum temperature for heating is 380°C. This model is widely used for heating in the school laboratories, industrial and mining, same as biochemistry, agriculture, petrochemicals engineering and pharmaceutical industry for the liquid.

Characteristics

- ✓ Shaping shell is made of cold rolled steel sheet with surface spraying treatment
- ✓ Nichrome heating wire, litter heating time, long service life
- ✓ Easy handling, wide application
- ✓ Precise PID temperature control meets the needs of high-precision experiments
- ✓ Digital type single position mantles have outside sensor, multi position mantles have inside sensor, optional outside sensor
- ✓ Series YR05021 to YR05029 are analogue, while the series YR05030 to YR05038 are type single position mantles with outer sensor, and the multiposition mantles are with inner sensor, and outer sensor optional.



#Let'sgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



KALSTEIN® Heating Mantles**Series YR05021**

Type	Analogical						
Model	YR05021	YR05022	YR05023	YR05024	YR05025	YR05026	YR05027
Capacity (ml)	100	250	500	1,000	2,000	3,000	5,000
voltage(v)	200 – 240V / 50 – 60Hz or 100 – 120V / 50 – 60Hz						
Maximum use temperature (°C)	380						
Power (W)	100	150	250	350	450	600	800
job size	Continuous						
Outer size (mm)	Ø200×160		Ø260×200		Ø300×230	Ø340×245	Ø350×250
Packing size (mm)	230x230x170		290x290x220		330x330x250	365x365x295	390x390x310
Net weight / kg)	2.5		4		5	7.5	8.5
Gross weight (kg)	2.8		6		7	8	9.2
Type	Digital						
Model	YR05030	YR05031	YR05032	YR05033	YR05034	YR05035	YR05036
Capacity (ml)	100	250	500	1,000	2,000	3,000	5,000
voltage(v)	200 – 240V / 50 – 60Hz or 100 _ 120V / 50 – 60Hz						
Maximum use temperature (°C)	380						
Power (W)	100	150	250	350	450	600	800
job size	Continuous						
Outer size (mm)	Ø200 x 160		Ø260×200		Ø300×230	Ø340×245	Ø350×250
Packing size (mm)	230x230x170		290x290x220		330x330x250	365x365x295	390x390x310
Net weight / kg)	2.5		4		5	7.5	8.5
Gross weight (kg)	2.8		6		7	8	9.2

**#Letsgivingmore**

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S





 **Kalstein**[®]
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



Fundación
Jacinto Convit

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



Humatem

Solidarité et coopération
biomédicale


<https://www.humatem.org/>



**MANIAPURE
FOUNDATION USA**

<https://www.maniapurefoundation.org/>

Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed

The 3 articles of the week!

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

1
Read this ARTICLE NOW



How does the climate affect your microscope?

microscope is an instrument that allows observing and measuring objects that are too small to be seen with the naked eye. This is achieved through the use of lenses or a system of mirrors that collect and focus light from the object being observed.

Kalstein Blog

Article of the week:
How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

2
Read this ARTICLE NOW



What is the importance of the laboratory refrigerator?

Refrigerator is a device that is one of the most common pieces of equipment that is used in laboratories. It is a controlled environment (refrigerated space), which allows and allows the storage of biological samples and reagents.

Kalstein Blog

Article of the week:
What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...

23 views • 3 years ago



Kalstein's Microscopes

445 views • 1 week ago



PARIS - FRANCE
HEADQUARTER

2 Rue Jean Lantier
75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89

Cel: +33 (0) 6 80 76 07 10

sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

