




# MEDICAL SECTOR

## INFANT PHOTOTHERAPY INCUBATOR



**Kalstein**®

A different accompaniment, at your service





# Kalstein®

A different accompaniment, at your service

OUR SERVICES

## Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



### Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



### Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



### #Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



### Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



### Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



### Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



**KALSTEIN® Infant Phototherapy Incubator**  
**Infant Phototherapy Incubator**  
**Model YR02184**



**#Letsgivemore** ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
2 Rue Jean Lantier • 75001 Paris •  
+33 1 78 95 87 89 / +33 6 80 76 07 10 •  
<https://kalstein.eu>  
**KALSTEIN FRANCE, S. A. S**



## KALSTEIN® Infant Phototherapy Incubator Infant Phototherapy Incubator

Model YR02184



### Features:

- ✓ Phototherapy light intensity is constantly adjustable
- ✓ Time of bilirubin phototherapy can be paused, recorded and accumulated
- ✓ Phototherapy treatment duration can be pre-set and it automatically shuts off when a pre-set duration is reached
- ✓ Working life of LED phototherapy unit is lasting 20,000 hours
- ✓ Silicone mattress protect neonate from came while phototherapy treatment
- ✓ Microprocessor servo-controller ensures a precise and reliable temperature control
- ✓ Set temperature, monitored air temperature and baby skin temperature and heating power are distinguished in display
- ✓ Servo-controlled humidification system maintains the relative humidity at a constant level in the incubator
- ✓ Set and monitored humidity concentration are displayed separately
- ✓ Ultrasonic and low temperature nebulized humidification protects neonate from bacteria's pollution
- ✓ >37°C override temperature function
- ✓ Self-checking function for safety affirmative
- ✓ Audible and visible alarm indications for multiple failures with 3 priorities
- ✓ Chord alarm notes nurse farther than 8 meters, but leaves asleep baby uninterrupted
- ✓ Aluminous framework for constant temperature keeping
- ✓ Continuous adjustment of mattress inclination up to 10°
- ✓ Drawer-type water reservoir is easily removable for cleaning and disinfection
- ✓ Nature air circulation system
- ✓ Cupboard cabinet
- ✓ RS-232 Connector

### Product Description:

As our top one infant incubator is not a simple connection of incubator and phototherapy unit. The sophisticated integration enable that newborn with hyper bilirubinirrubine could receive all around blue-ray phototherapy within an appropriate, sheltered and reliable environment. Every details including distance between phototherapy units and patient, illumination intensity and so on results from persistent and strict research and experience for your newborn patient while it is rigorously carrying out the duty of protecting infant with best care.

### Application:

The YR02184 not only can keep a stable environment for neonate patient, but also can treat baby jaundice with intensive LED phototherapy simultaneously. With its compact design, ergonomics height adjustment and smooth trolley, the caregivers are easy to use it friendly in NICU. Widely applicable in all hospitals and clinics and optimize for preterm births or for some ill full-term babies or jaundiced baby.

### Optional Accessories:

- ✓ Side door, front-rear door
- ✓ Double wall hood
- ✓ Infant phototherapy unit (OSRAM tubes)



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
2 Rue Jean Lantier • 75001 Paris •  
+33 1 78 95 87 89 / +33 6 80 76 07 10 •  
<https://kalstein.eu>  
KALSTEIN FRANCE, S. A. S



**KALSTEIN® Infant Phototherapy Incubator**  
**Infant Phototherapy Incubator**  
**Model YR02184**



**Technical Description:**

<b>Model</b>	YR02184
<b>Air temperature control</b>	•
<b>Control of skin temperature</b>	•
<b>37 ° air temperature control</b>	•
<b>The crib can be tilted</b>	•
<b>Humidity display</b>	•
<b>Humidity control</b>	•
<b>Superior Phototherapy Unit</b>	•
<b>Lower Phototherapy Unit</b>	•
<b>Side door</b>	○
<b>Aluminum tank for heat preservation.</b>	•
<b>Windows operation and access ports</b>	1 port, 4 windows
<b>Cabinet</b>	•
<b>Adjustable working height</b>	•
<b>Additional rear door</b>	○
<b>Baby Scale Bed</b>	○
<b>Silent DC motor</b>	○



**#Letsgivemore** ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
 2 Rue Jean Lantier • 75001 Paris •  
 +33 1 78 95 87 89 / +33 6 80 76 07 10 •  
<https://kalstein.eu>  
**KALSTEIN FRANCE, S. A. S**







 **Kalstein**<sup>®</sup>  
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE  
ACQUISITION  
OF A KALSTEIN  
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



**Fundación  
Jacinto Convit**

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



**Humatem**

Solidarité et coopération  
biomédicale

<https://www.humatem.org/>




**MANIAPURE  
FOUNDATION USA**

<https://www.maniapurefoundation.org/>



**Articles Kalstein**  
(Social Networks and Official Websites)



**Kalstein Blog**  
Not just be informed, be fully informed

**The 3 articles of the week!**

Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>



**1**  
Read this ARTICLE NOW

**How does the climate affect your microscope?**

**Kalstein Blog**

**Article of the week:**  
How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>



**2**  
Read this ARTICLE NOW

**What is the importance of the laboratory refrigerator?**

**Kalstein Blog**

**Article of the week:**  
What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>

**Kalstein Training Videos**  
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...

23 views • 3 years ago



Kalstein's Microscopes

445 views • 1 week ago



PARIS - FRANCE  
**HEADQUARTER**

2 Rue Jean Lantier  
75001 Paris - France

**Tlf:** +33 (0) 01 78 95 87 89

**Cel:** +33 (0) 6 80 76 07 10

[sales@kalstein.eu](mailto:sales@kalstein.eu)

<https://kalstein.eu/>  
<https://www.kalstein.fr/>

