



# LABORATORY SECTOR

REFRIGERATORS AND FREEZERS



**Kalstein**®

A different accompaniment, at your service





# Kalstein®

A different accompaniment, at your service

OUR SERVICES

## Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



### Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



### Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



### #Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



### Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



### Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



### Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



# KALSTEIN® Refrigerators and Freezers

## Low Temperature Freezer

Model YR05100



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
2 Rue Jean Lantier • 75001 Paris •  
+33 1 78 95 87 89 / +33 6 80 76 07 10 •  
<https://kalstein.eu>  
KALSTEIN FRANCE, S. A. S



# KALSTEIN® Refrigerators and Freezers

## Low Temperature Freezer

Model YR05100



### Product Description:

Covid Vaccine Refrigerators / Standing Deep Freezer 528 Liters with separated double door.

### Safety:

- ✓ Multiple malfunction alarms :high and low temperature , sensor error, power off, temperature controller malfunction, door ajar alarms.
- ✓ Alarm indications: audible buzzing alarm and visual flashing light alarm
- ✓ Corrosion proof cabinet interior design
- ✓ Power failure protect: the device would be attack by the high current when all devices restart simultaneously, the start-delay function would protect device from the damage

### Reliability:

- ✓ Chemically stable, CFC -free, commercially available and environmentally safe refrigerant.
- ✓ Optimized refrigeration system designed to maximize cooling effect and temperature uniformity
- ✓ Microprocessor control, digital display, adjustable temperature range: – 30°C–40°C
- ✓ LED digital display for clear observation,
- ✓ Wide voltage tolerance design with applicable voltage range of 187-242V/AC

### Ergonomic Design:

- ✓ USB data logging and temperature recorder (optional).
- ✓ Build-in lock with safety and beauty
- ✓ Lockable caster design for customer to move as needed

### Main Features:

- ✓ Microprocessor control, adjustable temperature range: -30°C to -40°C
- ✓ LED digital display and increment at 0.1°C
- ✓ Permanently lubricated cooling fan for safety and longevity .
- ✓ Independent inner door to avoid more cooling leaking out, and provide more stable inside temperature performance
- ✓ High density CFC-Free urethane foam insulation.



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
2 Rue Jean Lantier • 75001 Paris •  
+33 1 78 95 87 89 / +33 6 80 76 07 10 •  
<https://kalstein.eu>  
KALSTEIN FRANCE, S. A. S



**KALSTEIN® Refrigerators and Freezers****Low Temperature Freezer****Model YR05100****Technical Description:**

Power & Material		Accessories	
<b>Voltage/frequency(V/Hz)</b>	220/50Hz	Foot	Y
<b>Power(W)</b>	340W	Castor	Y
<b>Current(A)</b>	3.4A	Inside door	9 Drawer×2
<b>Inside material</b>	Color sprayed steel	Test hole QTY/diameter	2/25MM
<b>Outside material</b>	Color sprayed steel	Remote alarm port	Y
<b>Insulation</b>	PURF		
Dimension		Alarm System	
<b>Capacity(L/cu.ft)</b>	528/18.64	High/low temperature	Y
<b>NT./GT.(kg)</b>	151/184	Power failure	Y
<b>Exterior size(W*D*H)(mm)</b>	910*845*1875	Sensor error	Y
<b>Interior size(W*D*H)(mm)</b>	(651*631.5*635)*2	Temperature controller failure	Y
<b>Package size(W*D*H)(mm)</b>	930*970*2000	Power failure backup system(alarm)	8h
<b>CBM</b>	1.80	Certification	
<b>20GP/40GP/40HQ</b>	12/ 24/24	ISO9001,ISO13485,ISO14001,CE	

**#Letsgivemore** ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
2 Rue Jean Lantier • 75001 Paris •  
+33 1 78 95 87 89 / +33 6 80 76 07 10 •  
<https://kalstein.eu>  
**KALSTEIN FRANCE, S. A. S**





 **Kalstein**<sup>®</sup>  
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE  
ACQUISITION  
OF A KALSTEIN  
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



**Fundación  
Jacinto Convit**

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



**Humatem**

Solidarité et coopération  
biomédicale


<https://www.humatem.org/>



**MANIAPURE  
FOUNDATION USA**

<https://www.maniapurefoundation.org/>

**Articles Kalstein**  
(Social Networks and Official Websites)



**Kalstein Blog**  
Not just be informed, be fully informed

**The 3 articles of the week!**

Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>



**1**  
Read this ARTICLE NOW

**How does the climate affect your microscope?**

**Kalstein Blog**

**Article of the week:**  
How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>



**2**  
Read this ARTICLE NOW

**What is the importance of the laboratory refrigerator?**

**Kalstein Blog**

**Article of the week:**  
What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>

**Kalstein Training Videos**  
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...

23 views • 3 years ago



Kalstein's Microscopes

445 views • 1 week ago



PARIS - FRANCE  
**HEADQUARTER**

2 Rue Jean Lantier  
75001 Paris - France

**Tlf:** +33 (0) 01 78 95 87 89  
**Cel:** +33 (0) 6 80 76 07 10  
[sales@kalstein.eu](mailto:sales@kalstein.eu)

<https://kalstein.eu/>  
<https://www.kalstein.fr/>

