













All rights reserved ® KALSTEIN France S. A. S., 2 Rue Jean Lantier • 75001 Paris • +33 1 78 95 87 89 / +33 6 80 76 07 10 • https://kalstein.eu KALSTEIN FRANCE, S. A. S





OUR SERVICES Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings In any part of the world, receive your induction or training from our

induction or training from our specialized team of engineers



Quick Response Our work team is always available to response all your consults or questions, in order to support you in any situation.





#Letsgivemore Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or.air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.





KALSTEIN® Microscope Microscope Model YR05791-1

CE





#Letsgivemore 💗

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved ® KALSTEIN France S. A. S., 2 Rue Jean Lantier • 75001 Paris • +33 1 78 95 87 89 / +33 6 80 76 07 10 • https://kalstein.eu KALSTEIN FRANCE, S. A. S



All rights reserved © KALSTEIN France S. A. S.,

www.kalstein.eu

KALSTEIN® Microscope Microscope **Model YR05791-1**

Product Accessories



Abbe Condenser N.A.1.25 **Telescope 10x Green Filter**



Phase Contrast Objective 40x Phase Contrast Diaphragm 40x

Product Specification

	•
ltems	YR05791-1 Phase Contrast Microscope
Magnification	40x~1000x
Head	Sliding Binocular Head, 45° Inclined
Eyepiece	Wide Field WF10x/18mm
Nosepiece	Quadruple Nosepiece
Objective	Achromatic 4x/0.10, 10x/0.25, 100xs/1.25(0il)
Stage	Double Layer Mechanical Stage, 140*140mm, Mvoing Range 50*75mm
Focusing	Coaxial Coarse & Fine Focusing. Fine Focusing Scale 0.002mm, Coarse Focusing Stroke 37.7mm/Rotation, With Tension Adjustment, With UpLimit
Condenser	N.A.1.25 Abbe Condenser, Iris Diaphargm, With Filter Holder
Filter	Blue, Yellow, Green
Phase Contrast Illumination	Phase Contrast Objective 40x
	Phase Contrast Diaphragm 40x
	Abbe Condenser N.A.1.25
	Centering Telescope 10x
	Green Filter
	Plan-Concave Mirror
YR05791-1 Phase Contrast	Built-in Illumination 6V20W Halogen Lamp, 100~240V Microscope Optional Accessory
Head	Sliding Trinocular Head, 45° Inclined, Anti-Mould System
	Seidentopf Binocular Head, 30° Inclined, Anti-Mould System
	Seidentopf Trinocular Head, 30° Inclined, Anti-Mould System
Eyepiece	Wide Field Eyepiece P16x/11mm
Nosepiece	Quintuple Nosepiece



#Letsgivemore 💗 Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.





All rights reserved © KALSTEIN France S. A. S.,

www.kalstein.eu





#LETSGIVEMORE

WITH THE ADCQUISITION OF A KALSTEIN EQUIPMENT

YOU MAKE A CONTRIBUTION TO:





- C-- C-





ANY DOUBT?

Contact us!

Articles Kalstein (Social Networks and Official Websites)

& Kalstein Blog

Not just be informed, be fully informed

The 3 articles of the week!







Article of the week: How does the climate affect your microscope?



Kalstein Blog

ial medias or official websites https://kalstein.eu/articles

Article of the week: What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites https://kalstein.eu/articles



PARIS - FRANCE **H E A D Q U A R T E R** 2 Rue Jean Lantier 75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89 **Cel**: +33 (0) 6 80 76 07 10 sales@kalstein.eu

https://kalstein.eu/ https://www.kalstein.fr/

Kalstein Training Videos

(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...

23 views · 3 years ago



Kalstein's Microscopes445 views • 1 week ago

