




MEDICAL LINES

OPERATING TABLES



Kalstein®

A different accompaniment, at your service





Kalstein®

A different accompaniment, at your service

OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



Product overview

YR02154 is an electro-hydraulic operating table designed to meet universal surgical-positioning needs, like in Neurology, Ophthalmology, Gynecology, Urology, Minimal Invasive Surgery, Hand & Arm surgery.....

Its no-crossbar design allows for easy access to patient from both sides, 296-mm longitudinal shift provides bigger irradiated area but leaves minimal dead space, which makes YR02154 excellent solution for those procedures that incorporate X-ray imaging application, such as knee arthroscopy.

With corded remote or override control system, YR02154 can be set up smoothly and quickly to response to positioning variation during operation, or restored to the originally level position with just one press of reset button after surgical procedures end or before maintenance programs start.

To doubly ensure clinical safety, in addition to Auto-lock design of control system to avoid accidental mis-action, a rechargeable battery in YR02154 integrated and can be activated for backup of up to 50 batches of procedures, in case that AC power supply is cut off unexpectedly in operating rooms.

- Table body and sidebar are made from SUS-304 grade Stainless steel
- Comprehensive accessories option for different surgical positioning needs

Features



- Operation control with corded remote or override system, with Auto-lock of corded remote in 30s to avoid accidental mis-operation



- One press to restore originally "Level" position after surgery ends or for maintenance



- Manually controlled longitudinal shift up to 296 mm, gives easy access of device to patient



- Integrated X-Ray film cassette bracket under back section



- Doubly articulated Head plate to provide patient greater comfort in sleeping position



- Detachable Head section for connection of wrist support or head fixation (ophthalmology or neurology)
- Detachable Leg section for connection orthopedic extension



- Highly maneuverable trolley with front-ended brake to lock and release during movement



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
Optimum Business Center 450 Rue Baden Powell,
• 34000 Montpellier, France.
Tél: +33 467158849 / +33 680760710 / +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A. S



YR02154



Electro-hydraulic Operating Table

YR02154 is an electro-hydraulic operating table designed to meet universal surgical-positioning needs, like in Neurology, Ophthalmology, Gynecology, Urology, Minimal Invasive Surgery, Hand & Arm surgery.....



One-press reset to restore originally "Level" position.



Remote or override control, with Auto-lock in 30s to avoid accidental mis-operation.

YR02154 Specification

Table Length	2100 mm
Table Width	500 mm (not include side-rail)
Table Height (Raised/lowered)	1030 mm / 700 mm
Trendelenburg / Reverse-trendelenburg	30 ° / 30 °
Lateral Tilt (right/left)	20 ° / 20 °
Back Plate (up/down)	68 ° / 42 °
Head Plate (up/down)	57 ° / 72 °
Leg Plate (down/split)	96 ° / 65 ° x2
Flex / Reflex	220 ° / 115 °
Longitudinal Shift	296 mm
Patient Weight	200 kg (during positioning)
Input Power	AC 110/220V, 60/50Hz



- Detachable Head section for wrist support or head fixation (ophthalmology or neurology)
- Detachable Leg section for orthopedic extension



Longitudinal shift up to 296 mm, gives easy access of device to patient



Operating table tops are made of phenolic plastics. It can provide extremely high X-ray transmission



Highly maneuverable trolley with front-ended brake to lock and release during movement.



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S., Optimum Business Center 450 Rue Baden Powell, • 34000 Montpellier, France.
Tél: +33 467158849 / +33 680760710/ +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A. S





- Radio-lucent table top
- Replaceable table pad, latex-free and easy to clean.

YR02154 Specification

Table Length	2100 mm
Table Width	500 mm (not include side-rail)
Table Height (Raised/lowered)	1030 mm / 700 mm
Trendelenburg / Reverse-trendelenburg	30 ° / 30 °
Lateral tilt (right/left)	20 ° / 20 °
Back plate (up/down)	68 ° / 42 °
Head plate (up/down)	57 ° / 72 °
Leg plate (down/split)	96 ° / 65 ° x2
Flex / Reflex	220 ° / 115 °
Longitudinal shift	296 mm
Patient weight	200 kg (during positioning)
Input power	AC 110/220V, 60/50Hz

200KG



#Let'sgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
Optimum Business Center 450 Rue Baden Powell,
• 34000 Montpellier, France.
Tél: +33 467158849 / +33 680760710/ +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A. S



Position	Typical Accessories	Surgery Example
	 <p>Order No.: SSTFJ-60-00 - Lateral armboard (for prone, neuro, lateral and chair positions), padded - Mounting clamp included</p>	<ul style="list-style-type: none"> - Thoracic surgery - Kidney/Ureter surgery - Plastic surgery
	 <p>Order No.: SSTFJ-20-00 - Armboard with ball-joint, horizontal and vertical adjustment 360 degrees, padded, sold in pairs - Mounting clamps included</p>	
	 <p>Order No.: SSTFJ-40-00 - Lateral support (or called Body support, Lateral Brace), padded, sold in pairs. Mounting clamps included</p>	

Position	Typical Accessories	Surgery Example
	 <p>Order No.: SSTFJ-250-00 - Drain pan for gynecology and urology.</p>	<ul style="list-style-type: none"> - Gynecology - Urology - Rectal Surgery
	 <p>Order No.: SSTFJ-30-00 - Leg support (or called Knee crutch), padded, sold in pairs. - Mounting clamps included.</p>	

Position	Typical Accessories	Surgery Example
	 <p>Order No.: TJ10300-E - Optional Headrest (for neurological and ophthalmologic surgery)</p>	<ul style="list-style-type: none"> - Shoulder Surgery - Neurology - Breast reconstruction surgery
	 <p>Order No.: SSTFJ-50-00 - Shoulder support, padded, sold in pairs. - Mounting clamps included.</p>	

Position	Typical Accessories	Surgery Example
	 <p>Order No.: F-03 - Orthopedic extension device with 4-caster accessory cart (Traction height, distance and direction adjustable)</p>	<ul style="list-style-type: none"> - All surgeries involving posterior back (cervical spine, back, rectal area and dorsal extremities) - Knee arthroscopy
	 <p>Order No.: SSTFJ-100-00 - Easy armboard, horizontal adjustment 180°, padded.</p>	

Position	Typical Accessories	Surgery Example
	 <p>Order No.: SSTFJ-70-00 - Foot-rest, padded, sold in pairs. - Mounting clamps included.</p>	<ul style="list-style-type: none"> - Neck surgery - Laparoscopic surgeries in pelvic or lower abdominal region

Position	Typical Accessories	Surgery Example
	 <p>Order No.: SSTFJ-10-00 - L-shaped Anesthesia screen (or called Ether screen) with 2 wings, rigid - Mounting clamp included.</p>	<ul style="list-style-type: none"> - Spinal Surgery - Ophthalmology - Plastic Surgery - MIA (Minimal Invasive Access) surgery
	 <p>Order No.: SSTFJ-80-00 - External Kidney Elevator, adjustable height 0-120mm</p>	



#Letsgivemore ❤️

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
Optimum Business Center 450 Rue Baden Powell,
• 34000 Montpellier, France.
Tél: +33 467158849 / +33 680760710 / +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A. S



Gel pads and Positioner

For better access to surgical site and minimal risk of decubitus ulcer (or called "bedsore")

Prone Head Rest



Order No.	Description	Dimension	Quantity
HF-1101-1	Large	280x240x140 mm	1 pc
HF-1101-2	Medium	250x200x100 mm	1 pc
HF-1101-3	Small	230x200x115 mm	1 pc

Horseshoe-Shaped Head Rest



Order No.	Description	Dimension	Quantity
HF-1103-1	For Adult	200x70x70 mm	1 pc
HF-1103-3	For Child	140x50x30 mm	1 pc
HF-1103-5	For Neonate	80x30x20 mm	1 pc

Contoured Headrest



Order No.	Description	Dimension	Quantity
HF-1104	General purpose	180x180x80 mm	1 pc

Arm and Elbow pad



Order No.	Description	Dimension	Quantity
HF-2101-1	Medium	400x160x10 mm	1 pc
HF-2101-2	Large	600x160x10 mm	1 pc

Elbow protection strap



Order No.	Description	Dimension	Quantity
HF-2012	General purpose	400x150x15 mm	1 pc

Prone Chest support



Order No.	Description	Dimension	Quantity
HF-3102-1	PU material	630x460x150 mm	1 pc
HF-3102-3	Visco-elastic PU	600x490x180 mm	1 pc

Table cover pad (3 sections in one package)



Order No.	Description	Dimension	Quantity
HF-5101	Head section	500x230x15 mm	1 pc
	Back & Seat section	1000x500x15 mm	1 pc
	Leg section	500x500x15 mm	1 pc

Heel protection pad



Order No.	Description	Dimension	Quantity
HF-4103	Medium	150x80x45 mm	1 pc
HF-4103-1	Large	180x100x70 mm	1 pc

Ankle protection strap



Order No.	Description	Dimension	Quantity
HF-4105	General purpose	450x70x10 mm	1 pc

Ring-Shaped Head Rest



Order No.	Description	Dimension	Quantity
HF-1102-1	For Adult	200x70x70 mm	1 pc
HF-1102-3	For Child	140x50x30 mm	1 pc
HF-1102-5	For Neonate	80x30x20 mm	1 pc

Supine Head Rest



Order No.	Description	Dimension	Quantity
HF-1105-1	Large	230x230x40 mm	1 pc
HF-1105-2	Small	220x190x60 mm	1 pc
HF-1105-3	Medium	240x220x70 mm	1 pc

Headrest for Thyroid Gland Positioning



Order No.	Description	Dimension	Quantity
HF-1110-3	General purpose	480x390x90mm	1 pc

Lateral Body support



Order No.	Description	Dimension	Quantity
HF-3101-1	PU material	720x460x150 mm	1 pc
HF-3101-2	Visco-elastic PU	700x500x160 mm	1 pc

Wrist protection strap



Order No.	Description	Dimension	Quantity
HF-2013	General purpose	425x50x10 mm	1 pc

Contoured Armrest



Order No.	Description	Dimension	Quantity
HF-4102-1	Large	500x150x40 mm	1 pc
HF-4102-2	Medium	400x150x40 mm	1 pc

Rectangular Gel Pad for general purpose



Order No.	Description	Dimension	Quantity
HF-4104-1	Small	300x200x15 mm	1 pc
HF-4104-2	Medium	400x240x15 mm	1 pc
HF-4104-3	Large	500x300x15 mm	1 pc

D-shaped Gel pad for general purpose



Order No.	Description	Dimension	Quantity
HF-4101-1	Medium	300x150x70 mm	1 pc
HF-4101-3	Large	500x150x70 mm	1 pc
HF-4101-7	Filler for HF-4101-1		1 pc
HF-4101-8	Filler for HF-4101-3		1 pc



#Let'sgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
Optimum Business Center 450 Rue Baden Powell,
• 34000 Montpellier, France.
Titl: +33 467158849 / +33 680760710/ +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A. S



YR02155/6/7 series



Electro-hydraulic Operating Table

YR02155/6/7 Series operating tables are tailored to a variety of positioning requirements from different surgical disciplines. You always find the one from YR02155/6/7 Series that meet your demands perfectly, no matter during ENT, ophthalmological, endoscopic, or any other procedures.



YR02155/6/7 Series Specification

Model	YR02155	YR02156	YR02157
Control	Corded Remote Optional Pedal Control	Corded Remote N/A	Corded Remote N/A
Table Length	2100mm	2100mm	2100mm
Table Width	500mm(not include side rails)	500mm(not include side rails)	500mm(not include side rails)
Table Height	1040mm/690mm	1040mm/690mm	800mm/550mm
Trendelenburg/Reverse	25°/25°	25°/25°	25°/25°
Lateral Tilt (right/left)	20°/20°	20°/20°	20°/20°
Head Plate (up/down)	90°/90°	90°/90°	90°/90°
Back Plate (up/down)	80°/12°	80°/12°	80°/12°
Leg Plate(up/down/split)	10°/100°/100° x 2	10°/100°/100° x 2	10°/100°/100° x 2
Flex/Reflex	192°/100°	192°/100°	192°/100°
Kidney Elevator	0~120mm,build-in	0~120mm, external, optional	0~120mm, external, optional
Operation Load	300KG	200KG	200KG
Input Power	AC 110/220V, 60/50Hz	AC 110/220V, 60/50Hz	AC 110/220V, 60/50Hz



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
Optimum Business Center 450 Rue Baden Powell,
• 34000 Montpellier, France.
Tél: +33 467158849 / +33 680760710/ +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A. S




YR02155/6/7 Series- Design For Varying Requirements



Integrated kidney elevator 0~120mm height adjustment

For:

YR02155



Interchangeable leg and head section give bigger irradiated area for C-arm use.

For:

YR02155


YR02156



Optional pedal control for ENT, MIS and endoscopic procedures.

For:

YR02155



Extremely low table top about 55cm, give a comfortable height for ENT, neurosurgery and reduces strain when operating in a seated position.

For:

YR02157



Integrated backup battery optional, can support about 50 operations.

For:

YR02155

YR02156

YR02157



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
Optimum Business Center 450 Rue Baden Powell,
• 34000 Montpellier, France.
Tél: +33 467158849 / +33 680760710/ +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A. S



Select Guide



Functionality	Purpose & Application	Model
Longitudinal Shift	Adding more space for in C-arm X-ray imaging surgeries, such as cardiosurgery, orthopedic, ERCP, lithotomy.	YR02153, YR02154
Interchangeable Head and Leg Section	Interchangeable leg and head section give bigger irradiated area for C-arm use.	YR02153, YR02155, YR02156, YR02157
Low Height Position	Wide tabletop height adjustment range with low height position (550~600mm) is suitable in neurosurgery, ENT and plastic surgeries.	YR02153, YR02157
Double Back Sections	Kidney elevator structure can be reached by double back plates, used in kidney surgery.	YR02153
Kidney Elevator	Build-in or external kidney elevator make it easier for positioning in kidney surgery.	Build-in: YR02153, YR02155/6/7, YR02158 External: YR02154
Pedal Control	Providing more convenience and efficiency especially in laparoscopy and endoscopic surgery.	YR02153, YR02155
Override Control	Backup control panel on the column offers more safety in emergency.	YR02153, YR02154
Backlit Remote Handle	Making it easy to control the OT in dark OR ambient environment, such as during the endoscopic and laparoscopic surgeries.	YR02153
Motorized Back Plates	The back plates are electro-driven, suitable for lithotomy, neurosurgery, and other sitting position surgeries.	YR02153, YR02155
Valley-shaped Base	New design table base provides excellent stability and safety for surgeries, like neurosurgery, plastic surgery.	YR02153



#Letsgivemore


Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
Optimum Business Center 450 Rue Baden Powell,
• 34000 Montpellier, France.
Tél: +33 467158849 / +33 680760710 / +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A. S



Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed
The 3 articles of the week!
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



1
Read this ARTICLE NOW
How does the climate affect your microscope?
A microscope is an instrument that allows observing and measuring objects that are too small to be seen with the naked eye. This is achieved through a complex system of lenses and mirrors, which form an image of the object being observed.
Kalstein Blog
Article of the week:
How does the climate affect your microscope?
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



2
Read this ARTICLE NOW
What is the importance of the laboratory refrigerator?
A refrigerator in a laboratory is one of the most important pieces of equipment. It is used to maintain, in a controlled environment (refrigerated space), various fluids and substances.
Kalstein Blog
Article of the week:
What is the importance of the laboratory refrigerator?
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...

23 views • 3 years ago



Kalstein's Microscopes

445 views • 1 week ago



MONTPELLIER - FRANCE
HEADQUARTER
Optimum Business Center,
450 Rue Baden Powell,
34000 Montpellier - France.

Fax: +33 (0) 4 67 15 88 49
Tel: +33 (0) 6 80 76 07 10
sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

Official Websites for Latin America.

