




MEDICAL SECTOR

PATIENT AND FETAL MONITOR



Kalstein®

A different accompaniment, at your service





Kalstein®

A different accompaniment, at your service

OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



KALSTEIN® Patient and Fetal Monitor
Patients Monitors
Model YR02170 // YR02171 // YR02172



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



KALSTEIN® Patient and Fetal Monitor

Patients Monitors

Model YR02170 // YR02171 // YR02172



Product Description:

NICU Monitor. Devised specillay for neonatal. 8.4 inch touch screen,ECG (3 lead), Suntech NIBP, Masimo SpO2 , RESP, singal TEMP, HR/PR, Multiple display interface, respiratory asphyxia analysis, 112mm printer, lithium battery.

Features:

- ✓ 8.4 inch TFT touchscreen, light weight & compact design;
- ✓ Weak current detection up to 0.1mv of ECG ensures monitoring performance on neonates;
- ✓ Anti-movement technology using on SpO2 and NIBP ensures accurate reading;
- ✓ Unique adjustable respiration lead for raising neonatal patient monitoring accuracy;
- ✓ 112mm built-in printer maximum output 7 waveforms simultaneously



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
 2 Rue Jean Lantier • 75001 Paris •
 +33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



KALSTEIN® Patient and Fetal Monitor

Patients Monitors

Model YR02170 // YR02171 // YR02172



Technical Specifications:

Model	YR02170	YR02171	YR02172
ECG			
Input	3 leads or 5 leads selectable		
	3-conductor cable: RA; THE; LL or R: L; F		
	5-conductor cable: RA; THE; RL; LL, V or R; L; N; F; C		
Selection of potential clients	3 leads: 1,1, 1		
	5 leads: I, II, III, AVR, AVL, AVF, V		
ECG waveform	2 channels		
Win selection	x0.25.x0.5, x1, x2, Automatic		
Sweep speed	12mm / s, 25mm / s, 50mm / s		
Heart rate for newborns	15-350 BPM		
Heart rate for kids	15-100BPM		
Precision	1bpm or 1%, whichever is greater		
Resolution	1 bpm		
Filter	Diagnostic mode: 0.05-100Hz		
	Monitor mode: 0.5-40Hz		
	Surgical mode: 1.0-25Hz		
Protection	Supports voltage of 5000VAC / 50Hz in isolation		
	Against electrosurgical interference and defibrillation		
Alarm range	15-350 bpm		
ST segment detection	Measurement range: -2.0 mv to 2.0 mv		
	Alarm range: -2.0 mv to 2.0 mv		
Arrhythmia analysis	Yes		
Alarm	Yes, visual and audible alarm, recoverable alarm event		
Rhythm detection	Yes		
Apnea Analysis Diagram	display		
Calculation of drug dose	Yes		
Standby mode	Yes		
Breathing			
Method	Thoracic impedance		
Measuring range	Newborn: 0-150 rpm		
Apnea alarm	Yes, no alarm, 1 0s, 15 s, 20 s, selectable		
Alarm	Yes, audible and visual alarm, recoverable alarm events		
Method	Car manual		
Standby mode	Yes		
NIBP			



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
 2 Rue Jean Lantier • 75001 Paris •
 +33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



KALSTEIN® Patient and Fetal Monitor

Patients Monitors

Model YR02170 // YR02171 // YR02172



Method	Oscillometric		
Measurement modes	Manual / Automatic / Continuous		
Displayed parameters	Systolic, diastolic and MAP (mean arterial pressure)		
Automatic measuring time	Adjustable		
Unit of measurement	mmHg / kPa selectable		
Type of measure	Systolic, diastolic, medium		
Systolic pressure range	Newborn: 40-135 mmHg		
Diastolic pressure range	Newborn: 10-100 mmHg		
Medium pressure range	Newborn: 20-110 mmHg		
Precision	± 6 mmHg		
Resolution	1 mmHg		
Alarm	Systolic, diastolic, medium		
Overpressure protection	Double safety protection		
	SpO2		
Measuring range	0% -100%		
Resolution	0.01		
Precision	± 2% (70% -100%); unspecified (0-69%)		
Alarm range	0% -100%		
Pulse rate	Range: 20bpm-254bpm		
	Resolution: 1 bpm		
	Accuracy: ± 3 bpm		
	Alarm range: 0 ~ 254 bpm		
Additional characteristics	8.4 »touch screen, ECG, SpO2 (Massimo), RESP, NIBP (Suntech), TEMP & PR, & built-in thermal printer, 112mm	8.4 »touch screen, ECG, SpO2 (Massimo), RESP, NIBP, TEMP & PR	8.4 »touch screen, ECG, SpO2, RESP, NIBP, TEMP & PR



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
 2 Rue Jean Lantier • 75001 Paris •
 +33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S





 **Kalstein**[®]
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



**Fundación
Jacinto Convit**

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



Humatem

Solidarité et coopération
biomédicale


<https://www.humatem.org/>



**MANIAPURE
FOUNDATION USA**

<https://www.maniapurefoundation.org/>

Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed

The 3 articles of the week!

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



1
Read this ARTICLE NOW

How does the climate affect your microscope?

Kalstein Blog

Article of the week:
How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



2
Read this ARTICLE NOW

What is the importance of the laboratory refrigerator?

Kalstein Blog


Article of the week:
What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...
23 views • 3 years ago



Kalstein's Microscopes
445 views • 1 week ago



PARIS - FRANCE
HEADQUARTER
2 Rue Jean Lantier
75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89
Cel: +33 (0) 6 80 76 07 10
sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

