




# MEDICAL SECTOR

## PATIENT AND FETAL MONITORS



**Kalstein**®

A different accompaniment, at your service





# Kalstein®

A different accompaniment, at your service

OUR SERVICES

## Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



### Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



### Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



### #Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



### Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



### Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



### Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



**KALSTEIN® Patient and Fetal Monitor**  
**Patient Monitors**  
**Model YR02176**



**#Letsgivemore**

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
2 Rue Jean Lantier • 75001 Paris •  
+33 1 78 95 87 89 / +33 6 80 76 07 10 •  
<https://kalstein.eu>  
**KALSTEIN FRANCE, S. A. S**



## KALSTEIN® Patient and Fetal Monitor

### Patient Monitors

Model YR02176



#### Product Description:

- ✓ 8.4 inch touchscreen; Embedded handle; Able to use module; 1.6kg of weight; 5 hours of battery life; No-fan design; Multiple interfaces for different departments.
- ✓ 8.4 inch touchscreen;
- ✓ Embedded handle and 1.6kg of weight, easy to carry;
- ✓ Able to use module, anesthetic depth and EtCO2 monitor modules available;
- ✓ 5 hours of battery life, covers the requirement of transfer of patient;
- ✓ Vehicle power supply supported;
- ✓ No-fan design ensures low noise;
- ✓ Multiple interfaces for different departments.



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
 2 Rue Jean Lantier • 75001 Paris •  
 +33 1 78 95 87 89 / +33 6 80 76 07 10 •  
<https://kalstein.eu>  
 KALSTEIN FRANCE, S. A. S



## KALSTEIN® Patient and Fetal Monitor

### Patient Monitors

Model YR02176



#### Technical Description:

<b>Model</b>	YR02176
<b>Input</b>	3 leads or 5 leads selectable
	3-conductor cable: RA; LA LL or R: L; F
	5 Lead wire: RA: LA: RL: LL: V or R: LN: F: C
<b>Selection of potential clients</b>	3 leads: 1, 1, 111
	5 leads: 1, II, II, AR, AVL, AVF. V
<b>ECG waveform</b>	2 channels
<b>Win selection</b>	x0.25, x0.5.x1, x2.x4, Automatic
<b>Sweep speed</b>	12.5mm / s, 25mm / s, 50mm / s
<b>Heart rate</b>	for newborns 15-350BPM
<b>Heart rate</b>	for Adult/Pediatric 15-300BPM
<b>Precision</b>	1 bpm or 1%, whichever is greater
<b>Resolution</b>	1 bpm
<b>Filter diagnostic mode</b>	0.05-100 Hz
<b>Monitor mode</b>	0.5-40 Hz
<b>Surgical mode</b>	1.0-25Hz
<b>Protection</b>	Withstands 4000VAC / 50Hz voltage in isolation Against electrosurgical interference and defibrillation
<b>Alarm range</b>	Newborn: 15-350 BPM
	Adult / pediatric: 15-300 BPM
<b>ST segment detection</b>	Measurement range: -2.0 mv to 2.0 mv
<b>Alarm range</b>	-2.0m and 2.0mv
<b>Arrhythmia analysis</b>	Yes
<b>Alarm</b>	audible and visual alarm, recoverable alarm event
<b>Rhythm detection</b>	Yes
<b>Calculation of drug dose</b>	Yes
<b>Standby mode</b>	Yes
Breathing	



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
 2 Rue Jean Lantier • 75001 Paris •  
 +33 1 78 95 87 89 / +33 6 80 76 07 10 •  
<https://kalstein.eu>  
 KALSTEIN FRANCE, S. A. S



## KALSTEIN® Patient and Fetal Monitor

### Patient Monitors

Model YR02176



<b>Method</b>	Tooracl impedance
<b>Measuring range</b>	Newborn / Pediatric: 0-150 RPM
	Adult: 0-120 RPM
<b>Apnea alarm</b>	No alarm, 10s, 15s, 20s, 25s, 30s, 35s 40s selectable
<b>Alarm</b>	Audible and visual alarm, recoverable alarm events
NIBP	
<b>Method</b>	Oscillometric
<b>Measurement modes</b>	Manual / Automatic / Continuous
<b>Displayed parameters</b>	Systolic, diastolic and MAP (mean arterial pressure)
<b>Automatic measurement time</b>	Adjustable
<b>Unit of measurement</b>	mmHg / kPa selectable
<b>Type of measure</b>	Systolic, Diastolic. Piss
<b>Systolic pressure range</b>	Adult: 40-270 mmHg
	Pediatric: 40-200 mmHg
	Newborn: 40-135 mmHg
<b>Diastolic pressure range</b>	Adult: 10-210 mmHg
	Pediatric: 10-150 mmHg
	Newborn: 10-100 mmHg
<b>Medium pressure range</b>	Adult: 20-230 mmHg
	Pediatric: 20-165 mmH
	Newborn: 20-110 mmHg
<b>Precision</b>	± 5 mmHg
<b>Resolution</b>	1 mmHg
<b>Alarm</b>	Systolic, Diastolic, Mean
<b>Overpressure protection</b>	Double safety protection



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
 2 Rue Jean Lantier • 75001 Paris •  
 +33 1 78 95 87 89 / +33 6 80 76 07 10 •  
<https://kalstein.eu>  
 KALSTEIN FRANCE, S. A. S





 **Kalstein**<sup>®</sup>  
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE  
ACQUISITION  
OF A KALSTEIN  
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



**Fundación  
Jacinto Convit**

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



**Humatem**

Solidarité et coopération  
biomédicale


<https://www.humatem.org/>



**MANIAPURE  
FOUNDATION USA**

<https://www.maniapurefoundation.org/>

**Articles Kalstein**  
(Social Networks and Official Websites)



**Kalstein Blog**  
Not just be informed, be fully informed

**The 3 articles of the week!**

Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>



**1**  
Read this ARTICLE NOW

**How does the climate affect your microscope?**

**Kalstein Blog**

**Article of the week:**  
How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>



**2**  
Read this ARTICLE NOW

**What is the importance of the laboratory refrigerator?**

**Kalstein Blog**

**Article of the week:**  
What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>

**Kalstein Training Videos**  
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Alliances

23 views • 3 years ago



Kalstein's Microscopes

445 views • 1 week ago



PARIS - FRANCE  
**HEADQUARTER**

2 Rue Jean Lantier  
75001 Paris - France

**Tlf:** +33 (0) 01 78 95 87 89

**Cel:** +33 (0) 6 80 76 07 10

[sales@kalstein.eu](mailto:sales@kalstein.eu)

<https://kalstein.eu/>  
<https://www.kalstein.fr/>

