




LABORATORY LINES

POLARIZING MICROSCOPE



Kalstein®

A different accompaniment, at your service





Kalstein®

A different accompaniment, at your service

OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.





Polarizing Microscope Model YR0270

Kalstein's Polarizing Microscope YR0270.

With the features of easy operation, complete functions and economy, it is an ideal instrument in geology, petroleum, coal, chemicals and physics fields. It is also widely used in academic demonstration and scientific research.



Trinocular observing is available.



#Letsgivingmore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
Optimum Business Center 450 Rue Baden Powell,
• 34000 Montpellier, France.
Tél: +33 467158849 / +33 680760710 / +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A. S





Polarizing Microscope

Model YR0270

Kalstein's Polarizing Microscope YR0270.

With the features of easy operation, complete functions and economy, it is an ideal instrument in geology, petroleum, coal, chemicals and physics fields. It is also widely used in academic demonstration and scientific research.



Specifications

Viewing Head	Sliding Binocular Head at 45°, 360° Rotatable	•
	Sliding Trinocular Head at 45°, 360° Rotatable	
Eyepiece	Eyepiece WF10×/ 18 with Scale of Crosshair	•
	Eyepiece WF10×/ 18	•
Objective	Strain free Achromatic Objective 4×, 10×, 40×	•
Nosepiece	Toward Quadruple Nosepiece, Center Adjustable	•
Analyzer	Rotatable Analyzer with Gradation 0°-90°	•
Bertrand Lens	Bertrand Lens, Sliding in/ out of Optical Path	•
Optical Compensator	λ Slip (first class red)	•
	1/ 4 λ Slip	•
	Quartz wedge	•
Revolving Round Stage	Diameter Φ 160mm, Graduated in 1° Increments, Minimum Resolution 6' when using Vernier Scale	•
Condenser	Abbe Condenser with Iris Diaphragm & Filter	•
Focusing	Coaxial Coarse & Fine Adjustment, Range 28mm, Fine Division 0.002mm	•
Polarizer	Sliding in/ out of Optical Path, Located on the Top of Collector	•
Illumination	6V/ 20W Halogen Lamp, Brightness Adjustable	•

Note: • Standard Outfit, ◦ Optional



#Letsgivemore


Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
Optimum Business Center 450 Rue Baden Powell,
• 34000 Montpellier, France.
Tél: +33 467158849 / +33 680760710/ +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A. S



Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed
The 3 articles of the week!
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



1
Read this ARTICLE NOW



How does the climate affect your microscope?
A microscope is an instrument that allows observing and measuring objects that are too small to be seen with the naked eye. This is achieved through a complex system of lenses and mirrors, which form

Article of the week:
How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

2
Read this ARTICLE NOW



What is the importance of the laboratory refrigerator?
A refrigerator in a laboratory is one of the most important equipment for scientists to maintain, in a controlled environment (refrigerated space), various fluids and substances.

Article of the week:
What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...

23 views • 3 years ago



Kalstein's Microscopes

445 views • 1 week ago



MONTEPELLIER - FRANCE
SEDE

Optimum Business Center,
450 Rue Baden Powell,
34000 Montpellier - France.

Fax: +33 (0) 4 67 15 88 49

Tif: +33 (0) 6 80 76 07 10

sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

Official Websites for Latin America.

