


LABORATORY LINES AUTOCLAVES





All rights reserved ® KALSTEIN France S. A. S., Optimum Business Center 450 Rue Baden Powell, • 34000 Montpellier, France. TIF: +33 467158849 / +33 680760710/ +33 663810023 https://kalstein.eu KALSTEIN FRANCE, S. A. S



All rights reserved \circledast KALSTEIN France S. A. S.,



our services Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings In any part of the world, receive your

induction or training from our specialized team of engineers



Quick Response Our work team is always available to response all your consults or questions, in order to support you in any situation.





#Letsgivemore Thanks to your purchase, a donation will be made to a non-profit foundation that fights against

breast cancer and helps most

vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics


We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or.air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.





Portable Type Vertical Cylindrical Pressure SteamAutoclave Model YR05692



Application

YR series portable pressure steam autoclave is the equipment that disinfects and sterilizes theobjectswith saturated steam, its structure is vertical and heating way is electric heating. It is composedofsterilization room, shell, lid, steam generator, control system, display system and power-supplysystem. It is widely used in CSSD, lab, scientific research institution, pharmaceutical, food industryandissuitablefor items which are high temperature-resisting and high pressure-resisting including unwrappedor wrappedinstrument, dressing, rubber, liquid, plant-in instrument, porous item, hollow instrument to be sterilized.

Features

- 1. High quality stainless steel 304 material.
- 2. Coal or gas or electric heating system.Ç
- 3. Silicone seal.
- 4. Double scale indicator pressure gauge.
- 5. Double safety valve to ensure safe.
- 6. High pressure faucet drainage water.
- 7. Easy operate, safe and reliable.
- 8. Automatic cut off the power if the water is not enough.



#Letsgivemore 💗

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved (8) KALS TEIN France S. A. S., Optimum Business Center 450 Rue Baden Powell, • 34000 Montpellier, France. Tit: +33 467158849 / +33 680760710/ +33 663810023 https://aistein.eu KALSTEIN FRANCE, S.A. S



www.kalstein.eu



Safety Protection Device

- 1. Safe protection of water lacking.
- 2. Electric safety protective device.
- 3. Emergency exhaust valve.
- 4. Over temperature & over pressure auto-protection.
- 5. Over current protector.
- 6. Safety valve of auto release against over-pressure.

Main Configuration

- 1. Control: Microprocessor
- 2. Steam generator: Built-in steam generator, stainless steel 304
- 3. Chamber material: Stainless steel 304
- 4. Loading: Stainless steel drum.
- 5. Language: English
- 6. Voltage: 110/220V single phase, 50/60Hz



#Letsgivemore 💗

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved (8) KALSTEIN France S. A. S., Optimum Business Center 450 Rue Baden Powell, • 34000 Montpellier, France. Tit: +33 467158849 / +33 690760710/ +33 663810023 https://aistein.eu KALSTEIN FRANCE, S.A. S



www.kalstein.eu

Portable Type Vertical Cylindrical Pressure SteamAutoclave Model YR05692



Technical Parameter

Model	YR05692
Volume	18L
Chamber size	Ø300×330mm
Working temperature	126°C
Working pressure	0.165Mpa
Packing size	400×400×420mm
G.W/N.W	16kg/15kg
Power	AC220V/50Hz/2.0KW



#Letsgivemore 💜

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved (8) KALSTEIN France S. A. S., Optimum Business Center 450 Rue Baden Powell, • 34000 Montpellier, France. Tit: +33 467158849 / +33 680760710/ +33 663810023 https://aidein.eu KALSTEIN FRANCE, S.A. S

CE

www.kalstein.eu

All rights reserved © KALSTEIN France S. A. S.,





#LETSGIVEMORE

WITH THE **ADCQUISITION OF A KALSTEIN EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:





https://onetreeplanted.org

Humatem Solidarité et coopération biomédicale https://www.humatem.org/



https://www.maniapurefoundation.org/



ANY DOUBT? Contact us!







PARIS - FRANCE **H E A D Q U A R T E R** 2 Rue Jean Lantier 75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89 **Cel:** +33 (0) 6 80 76 07 10 sales@kalstein.eu


https://kalstein.eu/ https://www.kalstein.fr/

Kalstein Training Videos (Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...

23 views • 3 years ago



Kalstein's Microscopes 445 views • 1 week ago

