




LABORATORY LINES

SHAKER



Kalstein®

A different accompaniment, at your service





Kalstein®

A different accompaniment, at your service

OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.




KALSTEIN® Rotation Mixer YR05763



Features

- User friendly key operation interface, led timing time and speed display.
- A variety of centrifugal pipe clamps are optional, suitable for 1.5-50ml micro centrifugal pipe.
- Tilt angle of rotary table: 0-90 °, with scale display.
- It is equipped with timing function. The cultivation time can be set at will within the range of 0-99'59", the display screen displays the remaining time, and the sound alarm signal can be sent at the end of timing, or it can be set to operate continuously at irregular time.
- Long life motor drive, accurate speed, long life, maintenance free.
- It can be mixed in a mild way to avoid cell rupture and continuous mixing for a long time
- The easy to remove tube clamp allows easy and quick cleaning of spilled samples.
- The instrument can be placed in the refrigerator, freezer and other wet environment.
- A variety of centrifuge tube fixture is optional, applicable to 1.5-50ml micro tube.
- One key to start and stop, easy to operate, can run for a long time, smooth operation, low noise.
- It can be mixed in a gentle way to avoid cell rupture and continuous mixing for a long time at a uniform rate
- Easily removable tube clip allows simple and quick cleaning of spilled samples.
- The instrument can be run vertically up and down or rotated for different experiment requires.
- With digital timer and speed display(only 80 model)



Kalstein's YR05763 disk mixer is used for tissue culture, extraction and determination of sedimentation rate (including chemical determination of various blood samples), such as: blood sample mixing, radioimmunoassay, glucose, protein, DNA, cell, urine, enzyme and other wet chemical analysis. It is widely used in biochemical, biological, cell, immune, chemical, animal and plant analysis, clinical experiments and many other aspects. Its main feature is that they can be mixed in a mild way to avoid cell rupture and continuous mixing for a long time. The tilt angle of the turntable can be adjusted in the range of 0-90 °. It can also be placed in the incubator as required

Model	YR05763
Speed	10rpm ~ 80rpm
Tilt	0° ~ 90°
Timer	0~99'59"
Electric	AC 100 ~230V, 50/60Hz
Power	30W
Products size	W.204 x D.319 x H.300mm
N.W/ kg	4.5 kgs



#Letsgivemore ❤️

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S.





A different accompaniment, at your service

WITH THE ACQUISITION OF A KALSTEIN EQUIPMENT

YOU MAKE A CONTRIBUTION TO:

#LETSGIVEMORE 



Fundación Jacinto Convit

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



Humatem

Solidarité et coopération biomédicale

<https://www.humatem.org/>



MANIAPURE FOUNDATION USA

<https://www.maniapurefoundation.org/>

Thank you!

Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed
The 3 articles of the week!
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



1
Read this ARTICLE NOW
How does the climate affect your microscope?
A microscope is an instrument that allows observing and measuring objects at a very small scale. This is achieved through a complex system of lenses and mirrors, which form...

Kalstein Blog
Article of the week:
How does the climate affect your microscope?
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



2
Read this ARTICLE NOW
What is the importance of the laboratory refrigerator?
A refrigerator is a device that is one of the most common pieces of equipment used in laboratories. It is a controlled environment (refrigerated space), various fluids and substances...

Kalstein Blog
Article of the week:
What is the importance of the laboratory refrigerator?
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Alliances
23 views • 3 years ago



Kalstein's Microscopes
445 views • 1 week ago



PARIS - FRANCE
HEADQUARTER
2 Rue Jean Lantier
75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89
Cel: +33 (0) 6 80 76 07 10
sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

