




DENTAL SECTOR

DENTAL UNIT



Kalstein®

A different accompaniment, at your service





Kalstein®

A different accompaniment, at your service

OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



KALSTEIN® Dental Units
Strong Stability Dental Unit
Model YR06005



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



KALSTEIN® Dental Units
Strong Stability Dental Unit
Model YR06005



Multifunction Foot Control:

- ✓ Upgraded foot control ensures valves longer lifetime.
- ✓ Hand-free control bowl flush and clip filling avoid cross infection.

Multifunction Dental Assistant Unit:

- ✓ Flexible arm with three arises moves freely.
- ✓ High efficiency touch less control panel with 9 buttons.
- ✓ 3-way syringe with warm water.
- ✓ Air suction handle with regulator.
- ✓ Pre-position for built-in curing light.

Ergonomic Dental Chair Unit:

- ✓ Wide chair and U-shaped backrest for optimum patient comfort.
- ✓ Ergonomic seamless interior. Fine fiber leather is optional.
- ✓ Two armrests, swivel right arm rest for easy patient in and out.
- ✓ Double articulator headrest for children and disabled patients.

Comfortable Stool:

- ✓ Integrated back support and armrest good for relax of back, elbow and hands.
- ✓ Aluminum stool base makes ultra-stability.
- ✓ Nylon wheel moves quiet and smoothly.



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
 2 Rue Jean Lantier • 75001 Paris •
 +33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



KALSTEIN® Dental Units Strong Stability Dental Unit Model YR06005



Standard Accessories:

- ✓ 24V noiseless DC motor
- ✓ Sensor dental Led light
- ✓ Multifunctional assistant device
- ✓ Down instrument tray
- ✓ LED X-ray viewer
- ✓ 3 way syringe (cold / hot)
- ✓ High power suction (Air pressure)
- ✓ Saliva ejector (Water pressure)
- ✓ With main switch for water , air ,electricity
- ✓ Dentist stool RS03
- ✓ High quality air water tubes
- ✓ PU cushion, easy for clean
- ✓ Warm water system
- ✓ Pure water system
- ✓ Detachable and rotatable ceramic spittoon



Technical Description:

Model	YR06005
Voltage	100/110/120/220/230V, 50/60Hz
Air pressure	550-800 Kpa
Water pressure	200-400 Kpa
Net weight	170 Kg
Gross weight	200 Kg



#Letsgivemore ❤️

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S





#LETSGIVEMORE

**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**


YOU MAKE A CONTRIBUTION TO:



<https://onetreepanted.org>



Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed

The 3 articles of the week!

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



1
Read this ARTICLE NOW

How does the climate affect your microscope?

Kalstein Blog

Article of the week:
How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



2
Read this ARTICLE NOW

What is the importance of the laboratory refrigerator?

Kalstein Blog

Article of the week:
What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...

23 views • 3 years ago



Kalstein's Microscopes

445 views • 1 week ago



PARIS - FRANCE
HEADQUARTER

2 Rue Jean Lantier
75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89
Cel: +33 (0) 6 80 76 07 10
sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

