



LABORATORY LINES

TEST CHAMBERS



Kalstein®

A different accompaniment, at your service





Kalstein®

A different accompaniment, at your service

OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore ❤️

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.





Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.










ALLCOLD™ refrigeration technology



ALLSENS™ Microcomputer PID control technology



ALLFLOW™ Prefect air current cycling



#Letsgivingmore ♥
Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



TEST CHAMBERS




⊕ SUS304 stainless steel inner




⊕ Drainage design



⊕ PT 100 temperature sensor



⊕ Efficient isolation design



⊕ PRO-Bright™ Professional lighting system



⊕ Precision humidity sensor

⊕ Precision humidity sensor



#Let'sgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S





Operation chart:



1

Press the power switch



2

Press the set button | Set



3

Set time and temperature
(▶ Select shift key,
▼▲ and up and down key)




4

Set up to complete
Automatic operation of the
machine, arrive time will be
automatic alarm , and then
stop running.


Optional accessories:



1. Portable printer



2. LCD touch screen controller



3. ALLSENS™ software



4. Air interface



5. Powersupply interface



6. UV sterilization system



#Let'sgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



TEST CHAMBERS

Scope




Chinese medicine film



Pharmaceutical excipients



Drug injection



Blood products



IVD Reagents



Drug powder-injection



Drug GMP certification



#Letsgivingmore


Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



TEST CHAMBERS



#Let'sgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



Drug stability Test Chamber Strong Illumination Stability Test Chamber Comprehensive Drug Stability Test Chamber



Apply to pharmaceutical companies to the evaluation of drug failure, accelerated testing and long-term trials, pharmaceutical GMP certification, ideal for pharmaceutical stability testing equipment. have test hole and cap is made of silicone material in the right of the chamber, easy for product test.

According to Chinese pharmacopoeia requirements ,AS3100 of Australia , FDA of the United States, EN, ICH and other related standards. In accordance with pharmacopoeia requirements, must test the stability of drug in temperature, humidity ,lighting and so on integrated environment for limited time IQ, OQ, PQ and other related certification services. CE certification. Two-year warranty period.

Refrigeration Technology


World famous compressor:balanced refrigeration technology, auto defrost,multiple security system, long time running,environmental protection,high efficiency and energy saving.

Auto-defrost function: original rapid heat pump defrost technology, for continuous operation, effectively prevent frost because of continuous operation of evaporator.

Refrigeration effect manual control system: meet customer's terrible work environment.

Imported DuPont SUVA R134 a environmental refrigerant. For cleaner world, for better work condition.

To save energy, to reduce the cost.



#Letsgivemore ❤️

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S




Drug stability Test Chamber Strong Illumination Stability Test Chamber Comprehensive Drug Stability Test Chamber



Perfect Air Current Cycling

ALLFLOW™ air current cycling: perfect forced convection, maximum samples number of working room, minimum temperature recovery time after opening, world famous axial fan, perfect air current design, ensure an ideal result of samples culture.

ALLFLOW™ air current cycling: keep temperature continuous and stable, keep humidity continuous and stable, ensure perfect sample culture environment.




Humidification System

Original design of humidification system, prevent disadvantages (non uniformity, instability) of traditional atomizer design.

Automatic water intake, automatic water level control, ensure the precision of humidity control.

Stainless steel SUS304 internal surface and water tank, ensure the quality of atomization and uses persistent.



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



Drug stability Test Chamber

Strong Illumination Stability Test Chamber

Comprehensive Drug Stability Test Chamber



Programmable PID Control

- Integrated dot-matrix LCD display, Chinese and English subtitles, all the parameters display clearly at a glance. friendly design, display parameters: temperature, humidity, illumination, cycle, stage, running time, run/stop.
- Adaptive PID controller precisely control the temperature and humidity, prevent temperature soaring, keep working room temperature stable and uniform.
- User password protection, built-in multifunctional memory menu, connect to multiple devices (up to 16 units) by RS485 interface, real-time monitoring.
- Programmable multi-sections control, to meet customer needs, temperature, humidity and illumination control up to 30 sets. Time setting: 0-99h, 0-9999m, programs can be executed automatically cycles.
- Set operations with beep tips
- Remote-controlled operation, with programmable software (option).

Excellent Imported Temperature and Humidity Sensor

- Imported industrial temperature sensor PT100, short response time, low self-heating temperature.
- Imported humidity sensor, aerospace material no maintenance, precise humidity control level.



Optional Product Accessories:

1. Portable printer
2. LCD touch screen controller
3. Air interface
4. Power supply interface
5. ALLSENS™ software
6. UV sterilization system



#Letsgivingmore ♥


Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



Drug stability Test Chamber
Strong Illumination Stability Test Chamber
Comprehensive Drug Stability Test Chamber



#Letsgivemore ❤️
Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



Drug Stability Test Chamber Strong Illumination Stability Test Chamber Comprehensive Drug Stability Test Chamber



Product Model	Drug Stability Test Chamber			Strong Illumination Stability Test Chamber			Comprehensive Drug Stability Test Chamber		
	YR05350	YR05351	YR05352	YR05353	YR05353-1	YR05353-2	YR05354	YR05355	YR05356
Convection Mode	Forced Convection								
Control System	Thirty stages Microprocessor PID controller								
Temp. Range	0~65°C			10~65°C (with lighting)/0~60°C (without lighting)r					
Temp.resolution	0.1								
Temp. Fluctuation (10-40°C)	±0.5								
Temp. Uniformity (10-40°C)	±1	±1	±1.5	±1	±1	±1.5	±1	±1	±1.5
Humidity Range	Humidity Range: 50~90%RH, Humidity Fluctuation: ±3%RH						Humidity Range: 50~90%RH, Humidity Fluctuation: ±3%RH		
Illumination Range				0~6000LX (five stages adjustable)					
Working environment	Ambient temperature: 10~30°C, Humidity<70%								
Insulation materials	Imported environmental protection type material								
External Dimensions (H×W×D)	1410×650×680	1730×650×740	1700×745×930	1410×650×680	1730×650×740	1700×745×930	1410×650×680	1730×650×740	1700×745×930
Internal Dimensions (H×W×D)	760×510×390	1100×510×450	1050×600×640	760×510×390	1100×510×450	1050×600×640	760×510×390	1100×510×450	1050×600×640
Standard tray number	3	4	4	3	4	4	3	4	4
Interior Volume(L)	150	250	400	150	250	400	150	250	400
Interior materials	Sus304 stainless steel materials								
Standard Power (W)	1080	1100	1350	1080	1100	1350	1080	1100	1350
Power supply voltage	Single phase AC220V/50Hz								
Net Weight(KG)	107	135	158	107	145	175	107.5	145	175
Shipping Weight(KG)	132	162	186	132	172	203	133	172	203
Shipping Dimensions (H×W×D)	1610×750×830	1930×750×890	1900×840×1080	1610×750×830	1930×750×890	1900×840×1080	1610×750×830	1930×750×890	1900×840×1080



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S





 **Kalstein**[®]
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



**Fundación
Jacinto Convit**

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



Humatem

Solidarité et coopération
biomédicale

<https://www.humatem.org/>



**MANIAPURE
FOUNDATION USA**

<https://www.maniapurefoundation.org/>

Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed
The 3 articles of the week!
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



1
Read this ARTICLE NOW
How does the climate affect your microscope?
A microscope is an instrument that allows observing and measuring objects that are too small to be seen with the naked eye. This is achieved through a complex system of lenses and mirrors, which form...

Kalstein Blog
Article of the week:
How does the climate affect your microscope?
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



2
Read this ARTICLE NOW
What is the importance of the laboratory refrigerator?
A refrigerator is a device that is one of the most common pieces of equipment used in laboratories to maintain, in a controlled environment (refrigerated space), various fluids and substances...

Kalstein Blog
Article of the week:
What is the importance of the laboratory refrigerator?
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Alliances
23 views • 3 years ago



Kalstein's Microscopes
445 views • 1 week ago



PARIS - FRANCE
HEADQUARTER
2 Rue Jean Lantier
75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89
Cel: +33 (0) 6 80 76 07 10
sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

