




# LABORATORY SECTOR

VISCOMETER



## Kalstein®

A different accompaniment, at your service





# Kalstein®

A different accompaniment, at your service

OUR SERVICES

## Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



### Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



### Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



### #Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



### Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



### Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



### Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



**KALSTEIN® Viscometer**  
**Viscosity Measuring Analyzer Laboratory Viscometer**  
**Model YR05855 // YR05856**



**#Letsgivingore** ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
2 Rue Jean Lantier • 75001 Paris •  
+33 1 78 95 87 89 / +33 6 80 76 07 10 •  
<https://kalstein.eu>  
**KALSTEIN FRANCE, S. A. S**



## KALSTEIN® Viscometer Viscosity Measuring Analyzer Laboratory Viscometer Model YR05855 // YR05856



### Product Description:

- ✓ Test range from 1 to 2,000,000 mpa.s (cps)
- ✓ High accuracy
- ✓ Automatic detection
- ✓ Easy to operate

### Attachments:

- ✓ No. 1,2,3,4 Spindles
- ✓ No. 0 Spindle
- ✓ Adapter
- ✓ Temperature Sensor
- ✓ Protection Frame
- ✓ Foundation Plate
- ✓ Up-Down Frame
- ✓ English Manual
- ✓ Hexagonal Spanner
- ✓ Solid Wrench



### Technical Description:

Model	YR05855	YR05856
Measurement Range	1-2,000,000 mpa.s / cps	
Display	LCD	
Spindles Specification	No. 0, 1, 2, 3, 4	
	Note: The bigger the number of the Spindles, the bigger the measuring value	
Rotating Speed of Spindles	0.3 / 0.6 / 0.5 / 3 / 6 / 12 / 30 / 60 rpm	
Automatic Switch	Free selection of proper rotating speed or Spindles number.	
Stable Cursor	Strip-type cursor fully charged displayed when reading basically stable	
Measurement Tolerance	± 2%	
Power	AC 220V ± 10% 50 Hz ± 10%	
Working Condition	Temperature 5-35 degree	
	Relative humidity no more than 80%	
Outline Dimension (mm)	370×325×280 mm	
Packing (mm)	400×360×310 mm	
Net Weight (kg)	6.8 kg	
Gross Weight (kg)	9 kg	
Additional Information	0.3,0.6,1.5,3,6,12,30,60 rpm/min	
Mpa.S	1-2,000,000 cps	



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
2 Rue Jean Lantier • 75001 Paris •  
+33 1 78 95 87 89 / +33 6 80 76 07 10 •  
<https://kalstein.eu>  
KALSTEIN FRANCE, S. A. S





 **Kalstein**<sup>®</sup>  
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE  
ACQUISITION  
OF A KALSTEIN  
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



**Fundación  
Jacinto Convit**

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



**Humatem**

Solidarité et coopération  
biomédicale


<https://www.humatem.org/>



**MANIAPURE  
FOUNDATION USA**

<https://www.maniapurefoundation.org/>

**Articles Kalstein**  
(Social Networks and Official Websites)



**Kalstein Blog**  
Not just be informed, be fully informed  
**The 3 articles of the week!**  
Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>



**1**  
Read this ARTICLE NOW  
**How does the climate affect your microscope?**  
microscope is an instrument that uses a series of lenses to magnify images of small objects. This is a common tool in laboratories and research centers.  
**Kalstein Blog**  
**Article of the week:**  
How does the climate affect your microscope?  
Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>



**2**  
Read this ARTICLE NOW  
**What is the importance of the laboratory refrigerator?**  
Refrigerator is a device that is one of the most important equipment that is used to maintain, in a controlled environment (refrigerated space), various fluids and substances.  
**Kalstein Blog**  
**Article of the week:**  
What is the importance of the laboratory refrigerator?  
Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>

**Kalstein Training Videos**  
(Youtube Channel and Official Websites)



**Kalstein France - Jacinto Convit Foundation / Allianc...**  
23 views • 3 years ago



**Kalstein's Microscopes**  
445 views • 1 week ago



PARIS - FRANCE  
**HEADQUARTER**  
2 Rue Jean Lantier  
75001 Paris - France

**Tlf:** +33 (0) 01 78 95 87 89  
**Cel:** +33 (0) 6 80 76 07 10  
[sales@kalstein.eu](mailto:sales@kalstein.eu)

<https://kalstein.eu/>  
<https://www.kalstein.fr/>

