

LABORATORY SECTOR

SPECTROPHOTOMETER







All rights reserved ® KALSTEIN France S. A. S., 2 Rue Jean Lantier • 75001 Paris • +33 1 78 95 87 89 / +33 6 80 76 07 10 • https://kalstein.eu KALSTEIN FRANCE, S.A.S





OUR SERVICES Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings In any part of the world, receive your

induction or training from our specialized team of engineers



Quick Response Our work team is always available to response all your consults or questions, in order to support you in any situation.





#Letsgivemore 💙

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or.air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.





KALSTEIN® Spectrophotometer Visible Spectrophotometer Model YR01847





#Letsgivemore 💗

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved ® KALSTEIN France S. A. S., 2 Rue Jean Lantier • 75001 Paris • +33 1 78 95 87 89 / +33 6 80 76 07 10 • https://kalstein.eu KALSTEIN FRANCE, S. A. S



All rights reserved © KALSTEIN France S. A. S.,

www.kalstein.eu

KALSTEIN® Spectrophotometer Visible Spectrophotometer Model YR01847



Product Description:

The instrument uses a lead screw structure so that the instrument wavelength accuracy and wavelength resolution can be greatly improved.

Auto Setting Wavelength:

Users set wavelength automatically through arrow keys to avoid operation errors.

Numerical Keys:

With microprocessor controlled, all parameters of the instrument can be easily set by numerical keys.

8mm Thick Optical Base:

The instrument use a rigid die-cast aluminum base as its optical mount to ensure instrument stability and reliability.

Data Output

Equipped with USB port to connect with a PC to display spectrum scanning, kinetics and multi wavelength testing results on the screen through the optional software.

Standard Curve

Instruments can set up various standard curves according to customer's solutions and find the concentration of unknown solutions.

Model	YR01847		
Wavelength Range	320-1100nm		
Bandwidth	2nm		
Wavelength Accuracy	±0.5nm		
Wavelength Repeatability	≤0.2nm		
Wavelength Setting	Auto		
Photometric Accuracy	±0.2%T		
Photometric Repeatability	0.2%T		
Photometric Range	-0.3-3A, 0-200%T, 0-9999C		
Stability	± 0.001A/h @ 500nm		
Stray Light	≤0.05%T@220nm, 360nm		
Data Output Port	USB Port		
Printer Port	Parallel Port		
Display	128*64 Dots LCD		
Lamps	Tungsten Lamp		
Detector	Silicon Photodiode		
Power Requirements	AC 220V/50Hz or AC 110V/60Hz		
Dimension	460*360*225mm		
Weight	16kg		



#Letsgivemore 💗

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved ® KALSTEIN France S. A. S., 2 Rue Jean Lantier • 75001 Paris • +33 1 78 95 87 89 / +33 6 80 76 07 10 • https://kalstein.eu KALSTEIN FRANCE, S. A. S



www.kalstein.eu

KALSTEIN® Spectrophotometer Visible Spectrophotometer Model YR01847



ltem	Description	Quantity	Unit
1	Spectrophotometer	1	set
2	1cm Glass cuvette	4	pcs
3	Power cord	1	pcs
4	User's Manual	1	pcs
5	Dust Cover	1	pcs



#Letsgivemore Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.

All rights reserved © KALSTEIN France S. A. S.,



All rights reserved ® KALSTEIN France S. A. S., 2 Rue Jean Lantier • 75001 Paris • +33 1 78 95 87 89 / +33 6 80 76 07 10 • https://kalstein.eu KALSTEIN FRANCE, S. A. S



www.kalstein.eu





#LETSGIVEMORE

WITH THE ADCQUISITION OF A KALSTEIN EQUIPMENT

YOU MAKE A CONTRIBUTION TO:





https://onetreeplanted.org





ANY DOUBT? Contact us!



Articles Kalstein (Social Networks and Official Websites)





PARIS - FRANCE **HEADQUARTER** 2 Rue Jean Lantier 75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89 **Cel:** +33 (0) 6 80 76 07 10 sales@kalstein.eu


https://kalstein.eu/ https://www.kalstein.fr/

Kalstein Training Videos (Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...

23 views • 3 years ago



Kalstein's Microscopes 445 views • 1 week ago

