



# Shanghai Ro-Chain Medical Co.,Ltd

Add: Room 207, Building 1, Kangqiao Road 787, Pudong New District, Shanghai, 201315  
Tel: 0086-21-58055329 0086-13248293519 E-mail:seven@ro-chainmed.com Web: www.ro-chainmed.com

## Distributor Price List for Infant Care Products

Product Name	Model No.	Product Photo	Products Description & Technical Data	Unit Price (S/FOB)	Packing information
Infant Incubator	RC-BIN3000A		1. Air temp. servo-controlled by computer 2. Separate setting temp. and real temp. display windows 3. Various and self-check alarms 4. Removable humidity trough, easy to clean 5. With the transmission shelf and tray 6. RS-232 connector 7. Phototherapy (option)	\$868.00	108*64*112cm 7292kgs 6.7*6cm*92kgs Volume Weight is 154kgs
	RC-BIN3000AS		1. Air temp. and skin temp. servo-controlled by computer; 2. Separate setting temp. air temp. skin temp. display windows 3. Various and self-check alarms 4. With the transmission shelf and tray 5. RS-232 connector 6. Phototherapy (option)	\$1,186.00	108*64*112cm 7292kgs 6.7*6cm*92kgs Volume Weight is 154kgs
	RC-BIN3000B		1. Air temp. and skin temp. servo-controlled by computer; 2. Various and self-check alarms 3. Free step mattress tilting adjustment; 4. Removable humidity trough, easy to clean 5. >37°C temperature setting 6. With the transmission shelf and tray 7. Added with guardrail 8. RS-232 connector 9. With drawer and cabinet 10. Phototherapy (optional)	\$1,486.00	108*64*92cm/0.64cm 51kgs/65kgs 108*64*70cm/0.46cm 44kgs/55kgs Total :1.12cbm/120kgs Volume Weight is 224kgs
	RC-BIN3000BD		1. Air temp and skin temp servo-controlled by computer; 2. Various and self-check alarms 3. Free step mattress tilting adjustment; 4. Removable humidity trough, easy to clean 5. >37°C temperature setting 6. With the transmission shelf and tray 7. Added with guardrail 8. RS-232 connector 9. With drawer and cabinet 10. Double wall hood 11. Phototherapy (optional)	\$1,728.00	108*64*92cm/0.64cm 51kgs/65kgs 108*64*70cm/0.46cm 44kgs/55kgs Total :1.12cbm/120kgs Volume Weight is 224kgs
	RC-BIN3000B Plus		1. Air temp and skin temp servo-controlled by computer; 2. Various and self-check alarms 3. Free step mattress tilting adjustment; 4. Removable humidity trough, easy to clean 5. >37°C temperature setting 6. With the transmission shelf and tray 7. Added with guardrail 8. RS-232 connector 9. With drawer and cabinet 10. Has humidity display and control function 11. Phototherapy (optional)	\$2,068.00	108*64*92cm/0.64cm 51kgs/65kgs 108*64*84cm/0.56cm 44kgs/55kgs Total :1.12cbm/120kgs Volume Weight is 244kgs
	RC-BIN3000B Deluxe		1. Air temp and skin temp servo-controlled by computer; 2. Various and self-check alarms 3. Free step mattress tilting adjustment; 4. Removable humidity trough, easy to clean 5. >37°C temperature setting 6. With the transmission shelf and tray 7. Added with guardrail 8. RS-232 connector 9. With drawer and cabinet 10. Has humidity display and control function 11. LCD Display Screen 12. Phototherapy (optional)	\$2,368.00	108*64*92cm/0.64cm 51kgs/65kgs 108*64*84cm/0.56cm 44kgs/55kgs Total :1.12cbm/120kgs Volume Weight is 244kgs
	RC-BIN4000A		1. 7" inch LCD touch screen control panel; 2. Air temp and skin temp servo-controlled by computer; 3. Various and self-check alarms 4. Free step mattress tilting adjustment; 5. Removable humidity trough, easy to clean 6. >37°C temperature setting 7. Added with guardrail 8. RS-232 connector 9. With drawer and cabinet 10. Has humidity display and control function 11. Phototherapy is optional 12. Electronic lift is optional RC-BIN4000A: Based on above features, with additional X-ray cassette RC-BIN4000B: Based on BIN-4000A, with additional X2 concentration control RC-BIN4000C: Based on RC-BIN4000B, with scale weighing function	\$3,868.00	108*70*90cm, 0.64cm, 68.85kgs 108*70*94cm, 0.64cm 46.63kgs Total :1.32cbm/140kgs
	RC-BIN4000B		\$4,768.00	108*70*90cm, 0.68cm, 68.85kgs 108*70*94cm, 0.64cm 46.63kgs Total :1.32cbm/140kgs	
RC-BIN4000C	\$5,668.00		108*70*90cm, 0.68cm, 68.85kgs 108*70*94cm, 0.64cm 46.63kgs Total :1.32cbm/140kgs		
Transport Incubator	RC-BT100		Transport Incubator 1. Flat with low noise and long life-span 2. Temperature deviation and over temperature could be recorded in front panel directly 3. A second thermal cut-out function for more safety 4. Six opening windows, two side access operating windows can be covered 5. Double wall hood and with white LED lighting lamp 6. With two oxygen cylinder and oxygen flow meter 7. With normal trolley (ambulance trolley is optional)	\$4,986.00	108*70*98cm/110/65kgs 149*64*95cm/60/80kgs 1.12cbm/170kgs Volume Weight is 187.00kgs
Infant Radiant Warmer	RC-BRW3000A		1. Skin temp and manual control by micro-computer 2. Over temperature alarm 3. Available adjustment of radiator and mattress tilting 4. With one illumination light 5. APGAR timer 6. RS-232 connector	\$768.00	117*75*109cm 64kg/84kgs 0.97cbm/94kgs Volume Weight is 194.00kgs
	RC-BRW3000B		1. Skin temp and manual control by micro-computer 2. Separate setting temp. real temp. and timer display windows 3. Available adjustment of radiator and mattress tilting 4. Four self-check alarms 5. With one illumination light 6. RS-232 connector 7. With two sides fluorescent phototherapy	\$898.00	117*75*109cm 64kg/84kgs 0.97cbm/94kgs Volume Weight is 194.00kgs
	RC-BRW3000B New		1. Skin temp and manual control by micro-computer 2. Separate setting temp. real temp. and timer display windows 3. Available adjustment of radiator and mattress tilting 4. Four self-check alarms 5. With one illumination light 6. RS-232 connector 7. With two sides fluorescent phototherapy 8. With X-ray cassette	\$945.00	117*75*109cm 64kg/84kgs 0.97cbm/94kgs Volume Weight is 194.00kgs
	RC-BRW3000B Plus		1. Skin temp and manual control by micro-computer 2. Separate setting temp. real temp. and timer display windows 3. Available adjustment of radiator and mattress tilting 4. Four self-check alarms 5. With one illumination light 6. RS-232 connector 7. 2inch 8. Oxygen bath(6.3LX2)and Oxygen flow meter 9. Two sides fluorescent phototherapy is optional	\$1,485.00	117*75*110cm/0.24cm 6584kgs 0.97cbm/94kgs Volume Weight is 162.00kgs



## Shanghai Ro-Chain Medical Co.,Ltd

Add: Room 207, Building 1, Kangqiao Road 787, Pudong New District, Shanghai, 201315  
Tel: 0086-21-58055329 0086-13248293519 E-mail:severn@ro-chainmed.com Web: www.ro-chainmed.com

### Distributor Price List for Infant Care Products

Product Name	Model No.	Product Photo	Products Description & Technical Data	Unit Price (\$/FOB)	Packing information
	RC-BRW4000A		RC-BRW4000A features as below: 1.7 inch LCD touch screen control panel 2.Skin temp and manual control by micro-computer 3.Separate setting temp,real temp and timer display windows 4.Variable adjustment of radiator and mattress tilting 5.Four self-check alarms 6.With one illumination light 7.RS-232 connector 8.With X-ray cassette 9.LED phototherapy is optional 9. Electronic lift is optional BRW-4000B: Based on above features has Suction and Oxygen bottle/S.O2 stand Oxygen flow meter	\$2,186.00	148*86*130 122/155kgs 1.65cmm/155kgs Volume Weight is 250.00kgs
	RC-BRW4000B		\$2,486.00	148*86*130 122/155kgs 1.65cmm/155kgs Volume Weight is 250.00kgs	
Infant Phototherapy Unit	RC-BBP400A(4 tubes )		Mobile type, 4 tubes blue fluorescent blabs, Blab life time 1000h.	\$205.00	/
	RC-BBP500A(5 tubes)		Mobile type, 5 tubes blue fluorescent blabs, Blab life time 1000h.	\$308.00	/
	RC-B100 (LED)		<ul style="list-style-type: none"> <li>Use LED as light source;</li> <li>LED Digital Timer, display therapy time &amp; total time;</li> <li>Lifting &amp; Mobile Stand, radiant head rotating by angle of 30°/ 60°/ 90°.</li> </ul>	\$345.00	96*52*26cm 0.13cmm/20kgs
	RC-BBP1000B(LED)		Ceiling type, LED blabs, Blab life time 5000h.	\$615.00	59*45*24cm N.W: 4kgs G.W: 7kgs
	RC-BBP1000B(LED)R		Blue fluorescent phototherapy	\$205.00	
	RC-BBP2000B(LED)		1.LCD display indicates the accumulative running time and the current time 2.Wave length:400-500nm 3.Fit in with any infant bed and incubator 4.1700mw/cm²(10.5mw/cm²/min) 5. Phototherapy intensity can be adjustable 6. LED bulb life time is 5000h.	\$645.00	80.5*58*21cm 0.19cmm/20kgs
	RC-BBP3000B(LED)			\$625.00	
	RC-BBP4000B(LED)			\$1,468.00	
	RC-BL360D		1.Small size, light and high irradiation 2.LED Blue light as the radiant source is irradiating uniformly and efficient 3.Circles 360 double sides for more effective irradiation 4.Two phototherapy mode: up side phototherapy and down side phototherapy, can be use separately. 5.The electronic ballast makes the power supply adaptable and the power factor high. 6.Therapy time of phototherapy units from up and down side are displayed separately. 7.Automatically record each therapy time and total used therapy time. 8.Air temperature and skin temperature separately displayed. 9.Power off memory,alarm: key for alarm,key-locking and self-examination function. 10.Four alarm functions for safety: Power failure, probe failure, over temperature, fan failure.	\$2,485.00	120*70*96cm 0.81cmm/90kgs
Jaundice meter	RC-JD01		Transcutaneous jaundice meter/ jaundice detector/ Bilirubin meter, with battery	\$958.00	160x2x25cm/1.3kg
Digital baby scale	RC-05		1.Capacity:20kg/5g 2.Power:4*AA batteries 3.Product size:65x41x9cm	\$58.00	Carton size:70x46x52.5cm/5pcs N.W/G.W :15kg/12kg
	RC-CB551		1.Capacity:20kg/5g 2.Power:4*AA batteries 3.Product size:620*98*332mm	\$48.00	Carton size:660*360*355mm/3pcs G.W:7.5kg
	RC-EBSH		1.Capacity: 20kg/grad.-5g. 2.Unit: kg & lb 3.Use 4 x AA Batteries	\$75.00	Carton size:63*45*39cm/3pcs G.W:12.5kg
	RC-SCRY		1.Capacity: 20KG/10G 2.Weighing Unit:kg/lb/oz 3.Powered : 4pcs AA	\$50.00	Carton size:36*38*56cm/2pcs N.W/G.W :7.5KG/3.1KG



## Shanghai Ro-Chain Medical Co.,Ltd

Add: Room 207, Building 1, Kangqiao Road 787, Pudong New District, Shanghai, 201315  
Tel: 0086-21-58055329 0086-13248293519 E-mail:seven@ro-chainmed.com Web: www.ro-chainmed.com

### Distributor Price List for **Infant Care Products**

Product Name	Model No.	Product Photo	Products Description & Technical Data	Unit Price (\$/FOB)	Packing information
	<b>RC-EBSD</b>		1.Capacity: 20kg/5g 2.Units: kg/lb 3.Power: 4*AA batteries	<b>\$75.00</b>	Carton size 73*59*38cm/4pcs/28.5kg
	<b>RC-EBSL</b>		1.Capacity: 20kg/grad.5g 2.Units: kg/lb 3.Power supply: 4*AA batter	<b>\$75.00</b>	CTN Size: 66.7*54.1*34.4cm/3pcs G.W/N.W: 11.5kg/7.8kg
Infant incubator's Spare Part	<b>Drawer</b>			<b>\$45.00</b>	/
	<b>Cabinet</b>			<b>\$235.00</b>	/
	<b>Electric weighing function</b>			<b>\$295.00</b>	/
	<b>Electrical lift</b>			<b>\$375.00</b>	/
Infant warmer's Spare Part	<b>Phototherapy</b>			<b>\$185.00</b>	/
	<b>LED Phototherapy LED</b>			<b>\$545.00</b>	/
	<b>Suction</b>			<b>\$57.00</b>	/
	<b>X-ray tray</b>			<b>\$57.00</b>	/
	<b>Oxygen Bottle</b>			<b>\$94.00</b>	/
	<b>Oxygen Flow Meter</b>			<b>\$25.00</b>	/

**Remarks :**

1. All above products are of CE and ISO Certificate.
2. All above product price is quoted based on Exchange Rate between US and RMB is around 7. If above 5% vary, the price will vary accordingly.
3. All above product are quoted based on 1 unit per order. If more than 5 units per order, the unit price could be negotiable.
4. All above product's production lead-time is 7-10 working days after receiving payment.
5. All above product's quotation's valid date is 6 months upon quoted date.